# **RIDASCREEN® Nitrofurantoin (AMOZ)**

**Art. No. R3722**

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# RIDASCREEN<sup>®</sup> Nitrofuran (AMTZ)

## Brief information

RIDASCREEN<sup>®</sup> Nitrofuran (AMTZ) (Art. No.: R3722) is a competitive enzyme immunoassay for the quantitative determination of AMTZ in shrimp, fish, meat (bovine and porcine) and poultry (chicken and turkey).

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 96 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation: homogenization, derivatization, extraction, centrifugation, evaporation, defatting

Depending on the workflow in the laboratory, a 3-hour or an over-night derivatization may be performed. For highest precision, the over-night method is recommended. All specifications in this manual were determined using this sample preparation; details on the 3 hour method are published in the validation report.

Time requirement:

method I (short)	
sample preparation.....	3 h
derivatization .....	3 h
method II (over-night)	
sample preparation.....	3 h
derivatization .....	over-night (approx. 16 h)
test implementation .....	45 min

Limit of Detection: (corresponding to the standard substance)

shrimp .....	approx. 30 ng/kg
fish .....	approx. 40 ng/kg
meat (bovine) .....	approx. 40 ng/kg
meat (porcine) .....	approx. 65 ng/kg
poultry (chicken, turkey) .....	approx. 40 ng/kg

Detection capability (cc $\beta$ ): (corresponding to the standard substance)	shrimp .....	300 ng/kg
	fish .....	300 ng/kg
	meat (bovine) .....	300 ng/kg
	meat (porcine) .....	300 ng/kg
	poultry (chicken, turkey) .....	300 ng/kg
Recovery (corresponding to the standard substance; spiked samples)	shrimp .....	approx. 100 %
	fish .....	approx. 90 %
	meat (bovine) .....	approx. 110 %
	meat (porcine) .....	approx. 110 %
	poultry (chicken, turkey) .....	approx. 90 - 100 %

The specificity of the RIDASCREEN<sup>®</sup> Nitrofurantoin (AMOZ) test was determined by analyzing the cross-reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from the values determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the limit of detection and the recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

#### Specificity:

in buffer	NP-AMOZ (Standard substance).....	100 %
	Furaltadone .....	ca. 70 %
	Other Nitrofurantoin-antibiotics, derivatives, metabolites und parent substances (NP-AOZ, NP-AHD, NP-SEM, AOZ, AMOZ, AHD, SEM, Furazolidone, Nitrofurantoin, Nitrofurazone).....	< 1 %
	Chloramphenicol .....	< 1 %

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website [www.r-biopharm.com/products/food-feed-analysis](http://www.r-biopharm.com/products/food-feed-analysis).

## Related products

RIDASCREEN® Nitrofuran (AOZ) .....	(R3703)
RIDASCREEN® Nitrofuran (AHD).....	(R3713)
RIDASCREEN® Nitrofuran (SEM) .....	(R3715)
RIDA® Nitrofuran (AMTZ) Dotierlösung.....	(R3799)
RIDA® Nitrofuran (AOZ) Dotierlösung.....	(R3798)
RIDA® Nitrofuran (AHD) Dotierlösung .....	(R3796)
RIDA® Nitrofuran (SEM) Dotierlösung .....	(R3797)

### 1. Intended use

RIDASCREEN® Nitrofuran (AMTZ) is a competitive enzyme immunoassay for the quantitative determination of AMTZ in shrimp, fish, meat (porcine and bovine) and poultry (chicken and turkey).

### 2. General

Nitrofurans are synthetic broad-spectrum antibiotics, which are frequently used in animal production due to their excellent antibacterial and pharmacokinetic properties. They have also been used as growth promoters during the production of shrimps, poultry and pigs. Long term animal experiments have shown that the parent compounds and their metabolites have carcinogenic and mutagenic characteristics. This led to the prohibition of nitrofurans for the treatment of animals used for food production. In 1993, the EU banned the nitrofurans furaltadone, nitrofurantoin and nitrofurazone for use in animals used as sources of food, and in 1995 the use of furazolidone was also prohibited. The analysis of nitrofurans are based on the detection of the tissue bound metabolites of nitrofurans. The parent compounds are difficult to detect accurately since they are metabolized very rapidly after treatment. The tissue bound nitrofuran metabolites however are present for a long time after administration and they are used to detect nitrofuran abuse.

Parent compound	Metabolite	Abbreviation
Nitrofurantoin	1-Aminohydantoin	AHD
Furaltadone	3-Amino-5-morpholinomethyl-2-oxazolidinone	AMTZ
Furazolidone	3-Amino-2-oxazolidinone	AOZ
Nitrofurazone	Semicarbazide	SEM

Prior to analysis, the metabolites have to be derivatized by incubation with 2-Nitrobenzaldehyde into NP-AHD, NP-AMTZ, NP-AOZ and NP-SEM.

### 3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies directed against anti-AMOZ antibodies. The microtiter plate is already functionalized with anti-AMOZ antibodies, which are bound by the immobilized capture antibodies. Standards or sample and AMOZ conjugate are added into the wells. Free and enzyme conjugated AMOZ compete for the antibody binding sites (competitive enzyme immunoassay). Any unbound conjugate is then removed in a washing step.

The conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorption is inversely proportional to the AMOZ concentration in the sample.

## 4. Reagents provided

Each kit contains sufficient materials for 96 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format	Volume
Microtiter plate	-	Ready-to-use	96 wells
Standard 1	White	Ready-to-use	0 ng/L 1,3 ml
Standard 2	White	Ready-to-use	125 ng/L 1,3 ml
Standard 3	White	Ready-to-use	300 ng/L 1,3 ml
Standard 4	White	Ready-to-use	750 ng/L 1,3 ml
Standard 5	White	Ready-to-use	1500 ng/L 1,3 ml
Standard 6	White	Ready-to-use	3000 ng/L 1,3 ml
Wash buffer salt Tween		Salt for dissolving	
Conjugate	Red	Ready-to-use	7,5 ml
Substrate/Chromogen Red Chromogen Pro	Brown	Ready-to-use	13 ml
Stop solution	Yellow	Ready-to-use	14 ml
2-Nitrobenzaldehyde		Solid	100 mg

## 5. Materials required but not provided

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For the derivatization step, a solution of 10 mM 2-Nitrobenzaldehyde (see 4.) in Dimethylsulfoxide (DMSO) has to be prepared directly before use (i.e. dissolve 7.6 mg 2-Nitrobenzaldehyde in 5 ml DMSO).

## 6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

## 7. Storage instructions

Store the kit at 2 - 8 °C. Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C.

The substrate/chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after expiry of the kit (see kit label).

Do not interchange individual reagents between kits of different lot numbers.

## 8. Indication of deterioration of reagents

- any bluish coloration of the reddish substrate/chromogen prior to test implementation
- a value of less than 0.6 absorbance units ( $A_{450\text{ nm}} < 0.6$ ) for the zero standard

## 9. Preparation of Samples

The samples should be stored in a cool place, protected against light.

### 9.1. Shrimps (3h and over-night sample preparation)

- homogenize the sample, use stomacher or mixer
- prepare 10 mM 2-Nitrobenzaldehyde in Dimethylsulfoxide (DMSO) (see 5.2.)
- mix 1 g of the homogenized sample with 4 ml distilled water, 0.5 ml 1 M HCl and **100**  $\mu\text{l}$  10 mM 2-Nitrobenaldehyde in DMSO (see 5.2) by shaking thoroughly

#### Derivatisation

- incubate for 3h at 50 °C (122 °F) or over-night (approx. 16 h) at 37 °C (98.6 °F)
- add 5 ml 0.1 M  $\text{K}_2\text{HPO}_4$ , 0.4 ml 1 M NaOH and 10 ml ethyl acetate, shake vigorously for 30 sec
- centrifuge: 10 min / 4000 g / at room temperature (20 - 25 °C / 68 - 77 °F)
- transfer 5 ml of the ethyl acetate layer (upper layer) into a new vial and evaporate to dryness at 60 °C (140 °F)
- dissolve the residue in 1 ml n-hexane and mix properly with 1 ml wash buffer for 30 sec
- centrifuge: 10 min / 4000 g / at room temperature (20 - 25 °C / 68 - 77 °F)
- use 50  $\mu\text{l}$  of the lower, aqueous phase per well in the assay

### 9.2. Fish (3h and over-night sample preparation)

- homogenize the sample, use stomacher or mixer
- prepare 10 mM 2-Nitrobenzaldehyde in Dimethylsulfoxide (DMSO) (see 5.2.)
- mix 1 g of the homogenized sample with 4 ml distilled water, 0.5 ml 1 M HCl and **200**  $\mu\text{l}$  10 mM 2-Nitrobenaldehyde in DMSO (see 5.2) by shaking thoroughly

## Derivatisation

- incubate for 3h at 50 °C (122 °F) or over-night (approx. 16 h) at 37 °C (98.6 °F)
- add 5 ml 0.1 M K<sub>2</sub>HPO<sub>4</sub>, 0.4 ml 1 M NaOH and 10 ml ethyl acetate, shake vigorously for 30 sec
- centrifuge: 10 min / 4000 g / at room temperature (20 - 25 °C / 68 - 77 °F)
- transfer 5 ml of the ethyl acetate layer (upper layer) into a new vial and evaporate to dryness at 60 °C (140 °F)
  - if the sample has a high fat content, fatty residues that cannot be evaporated may occur. In this case the evaporation step is completed and sample preparation can be continued
- dissolve the residue in 1 ml n-hexane and mix properly with 1 ml wash buffer for 30 sec
- centrifuge: 10 min / 4000 g / at room temperature (20 - 25 °C / 68 - 77 °F)
- use 50 µl of the lower, aqueous phase per well in the assay

### 9.3. Meat and poultry (3h and over-night sample preparation)

- homogenize the sample, use stomacher or mixer
- prepare 10 mM 2-Nitrobenzaldehyde in Dimethylsulfoxide (DMSO) (see 5.2.)
- mix 1 g of the homogenized sample with 4 ml distilled water, 0.5 ml 1 M HCl and **300** µl 10 mM 2-Nitrobenzaldehyde in DMSO (see 5.2) by shaking thoroughly

## Derivatisation

- incubate for 3h at 50 °C (122 °F) or over-night (approx. 16 h) at 37 °C (98.6 °F)
- add 5 ml 0.1 M K<sub>2</sub>HPO<sub>4</sub>, 0.4 ml 1 M NaOH and 10 ml ethyl acetate, shake vigorously for 30 sec
- centrifuge: 10 min / 4000 g / at room temperature (20 - 25 °C / 68 - 77 °F)
- transfer 5 ml of the ethyl acetate layer (upper layer) into a new vial and evaporate to dryness at 60 °C (140 °F)
- dissolve the residue in 1 ml n-hexane and mix properly with 1 ml wash buffer for 30 sec
- centrifuge: 10 min / 4000 g / at room temperature (20 - 25 °C / 68 - 77 °F)
- use 50 µl of the lower, aqueous phase per well in the assay

## 10. Test implementation

### 10.1. Preliminary comments

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

As **wash buffer** a PBS tween buffer is needed. Please use the wash buffer salt (see 4.) contained in the kit. Dissolve the entire buffer salt in one liter of distilled water. The ready to use washing buffer expires after approx. 4 - 6 weeks at 2 - 8 °C (36 - 46 °F).

Alternative: Dissolve the contents of the pouch in 100 ml of distilled water to obtain a 10fold concentrated washing buffer. This solution expires after approx. 8 - 12 weeks, store at room temperature (20 - 25 °C). Use 1 part 10fold concentrate and dissolve with 9 parts of distilled water to obtain the ready to use washing buffer.

### 10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 50 µl of each standard or prepared sample to separate duplicate wells.
3. Add 50 µl of the conjugate to the bottom of each well, mix gently by shaking the plate manually and incubate for 30 min at room temperature (20 - 25 C) in the dark.
4. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl washing buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
5. Add 100 µl of substrate/chromogen solution to each well. Mix gently by shaking the plate manually and incubate for 15 min at room temperature (20 - 25 °C) in the dark.
6. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 30 minutes after addition of stop solution.

## 11. Results

A special software, the RIDA<sup>®</sup>SOFT Win.NET (Art. No. Z9996), is available to evaluate the RIDASCREEN<sup>®</sup> enzyme immunoassays.

The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = \% \text{ absorbance}$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates on semilogarithmic graph paper against the AMOZ concentration [ng/L]. Standard concentrations were adjusted to compensate for matrix effects.

In order to obtain the AMOZ concentration in ng/L / ng/kg (ppt) actually contained in a sample, the concentration read from the calibration curve must be further multiplied by the corresponding dilution factor. When working in accordance with the regulation stated, the **dilution factor is 2**.

**For further information or applications please contact your local distributor or R-Biopharm AG.**

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

### R-Biopharm AG

Postanschrift / Postal Address:

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64297 Darmstadt, Germany

Sitz / Corporate Seat: Pfungstadt

Tel.: +49 (0) 61 51 - 81 02-0

Fax: +49 (0) 61 51 - 81 02-40

E-mail: [info@r-biopharm.de](mailto:info@r-biopharm.de)

[www.r-biopharm.com](http://www.r-biopharm.com)

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are registered trademarks of R-Biopharm AG