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<b>Beispiel: Konelab / Arena</b>			
<b>Test Definition:</b>			
Full name	Enzytec Iron		
Online name	#		
Test used	yes		
Test type	Photometric		
Result unit	mg/l		
Number of decimals	1		
Temperature	37 °C		
Test limit	Low 0	High 20	Unit mg/l
Initial absorbance	0.000	2.000	E
Dilution limit	0	20	mg/l
Secondary dilution 1 +	#	#	
Acceptance	Automatic		
Dilution 1 +	0.0		
Sample type	Wine (Other)		
Correction factor	1.00		
Correction bias	0.00		
<b>Calibration parameter</b>			
Calibration type	Linear		
Curve direction	Ascending		
Repeat time (d)	1	Err. (mA)	#
Points/calibrator	simple	Err. (%)	#
Acceptance	Automatic		
Response limit min/max	# / #		
Bias correction in use	No		
Type of calibrators	Series (separate is also possible)		
Calibrator identification	Iron calibrator (20 mg/l)		
Concentrations	0 to 10 mg/l		
<b>Distribution</b>			
Additional blank	YES (or fixed time)	Antigen excess	No
Reagent 1	ENZYT Iron R1		
volume (µl)	150		
dispense with	Water	Volume (µl)	10
wash reagent	None		
Sample volume (µl)	15	(adjust to sensitivity needs)	
dispense with	Water	Volume (µl)	10
wash reagent	None		
Incubation (sec)	180		
Measurement BLANK	End-point		
response min./max	# / #		
Reagent 2	ENZYT Iron R2		
volume (µl)	25		
wash reagent	None		
dispense with	Water	Volume (µl)	10
Incubation (sec)	300		
Measurement	End-point		
λ 1 / λ 2 (nm)	600 / 700	575 / 700 is also possible	
measurement type	Fixed timing		

<b>Beispiel: Lisa 200</b>		
Full Test name	Enzytec Iron Ferene	
Short name	Iron	
Units	mg/l	
Type of test	End-point calibration	
Filter (nm)	580	
1st reading = zero	NO (A1 must be measured for each sample)	
Waiting time 1 (tours before addition of R2)	4 (4 tours = 4 x 24 sec.)	
Number of measurements (= incubation time)	15 (15 tours x 24 sec = 6 min)	
Reagent 1	VOL (µl)	250
	DIL (µl)	0
	POS	#
Reagent 2	VOL (µl)	45
	DIL (µl)	0
	POS	#
Sample.	VOL (µl)	20 - 50 µl (adjust to sensitivity needs)
	DIL	0
Starter (of reaction)	REAGENT 2	
Waiting time 2 (tours after addition of R2)	0 (measurements start immediately after addition of R2)	
Calibration	1 degree (= linear calibration; 2 degrees = non-linear calibration is also possible)	
Blank = Standard	YES (water used for blank is used as 1 <sup>st</sup> point of the calibration curve)	
Number of standards	3	
Standard 1	VAL (mg/l)	2 (= example with 2 mg/l)
	POS	# (position on the instrument)
Standard 2	VAL (mg/l)	5 (= example with 5 mg/l)
	POS	#
Standard 3	VAL (mg/l)	10 (= example with 10 mg/l)
	POS	#
Nb repeat Sample/Control	1 (number of repetitions ; it is possible to choose 2)	
Control	VAL (mg/l)	# (concentration of QC sample)
	POS	# (position on the instrument)
	DEV	# (deviation accepted)
Pre-dilution	1	
Post-dilution	4	
Sample diluents	Physiologic solution	
Type rinsing	3	
Normal values	UPPER	20 (in mg/l)
	LOWER	0
Limit linearity	20 (in mg/l)	
Blanc OD limit	LOWER	0
	UPPER	500

# value entered by the operator

<b>Beispiel: Hitachi 717</b>	
TEST	[IRON]
ASSAY CODE	[2POINTS] - [24] - [50]
SAMPLE VOLUME (µl)	[30] - [5] (adjust to sensitivity needs)
REAGENT VOLUME R1	[250] - [50] - [NO]
REAGENT VOLUME R2	[40] - [20] - [NO] (can be lower)
WAVELENGTH	[700] [600]
CALIBRATION TYPE	[LIN] [1] [4]
STD.(1) CONC. POS.	[0] - [1] (example)
STD.(2) CONC. POS.	[2] - [2] (2 mg/l =example)
STD.(3) CONC. POS.	[5] - [3] (5 mg/l =example)
STD.(4) CONC. POS.	[10] - [4] (10 mg/l =example)
SD LIMIT	[0.1]
DUPLICATE LIMIT	[#]
SENSITIVITY LIMIT	0
LIMITE D.O (CINET)	[0] - [INCR.]
PROZONE CHECK	[32000] - [1]
EXPECTED VALUE	[#] - [#]
PANIC VALUE	[#] - [#]
INSTRUMENT FACTOR	[1.0]