RIDA[®] Ochratoxin A column

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RIDA[®] Ochratoxin A column

Brief information

RIDA[®] Ochratoxin A column (Art. No.: R1303) are immunoaffinity columns for sample clean up prior to analysis of Ochratoxin A in food and feed.

R1303 contains 10 immunoaffinity columns for single use.

Gel suspension:	Sepharose
Antibody:	monoclonal antibody conjugated to Sepharose
Max. applied volume:	60 ml
Flow rate:	1 drop/sec
Detection limit:	depending on sample volume applied and on detection method
Capacity of the column:	approx. 200 ng ochratoxin A

1. General

RIDA® Ochratoxin A column are usable in combination with the enzyme immunoassays (R1311, RIDASCREEN® Ochratoxin A 30/15 and R5402, RIDASCREEN®FAST Ochratoxin A) for the quantitative determination of Ochratoxin A.

The sample preparation with the immunoaffinity columns simplifies and enhances the sample clean up procedure. Pure extracts are obtained, which can be analyzed by different analytical methods.

2. Test principle

The basis is the antigen-antibody reaction.

The column contains a gel suspension to which monoclonal antibodies are attached covalently. The antibodies are specific for the Ochratoxin A.

- A: The sample is applied and is passed through the column. If Ochratoxin A is present in the sample, it will be bound by the monoclonal antibodies.
- B: All other substances will not be retained on the column.
- C: Using methanol as eluent the Ochratoxin A will be released from the antigen-antibody-complex. Methanol causes a denaturation of the antibodies. Therefore, the antigen (Ochratoxin A) is set free and can be eluted.


Fig. 1: Principle of an immunoaffinity column

3. Reagents provided

With the immunoaffinity columns 10 sample preparations respectively could be done (one sample preparation per column). Each test kit contains:

10 immunoaffinity columns

4. Material required but not provided

4.1. Equipment

- laboratory balance
- grinder (mill) or coffee mill
- centrifuge
- shaker
- magnetic stirrer
- filter paper
- one way plastic syringes (as sample reservoir)
- optional: vacuum unit for simultaneous use of several columns
- graduated pipettes
- pasteur pipettes
- variable 20 µl - 200 µl and 200 - 1000 µl micropipettes

4.2. Reagents

- methanol
- acetonitrile
- n-heptane (or hexane)
- 20 mM PBS/methanol, pH 7.4 (90/10 (v/v)):
0.55 g $\text{NaH}_2\text{PO}_4 \times \text{H}_2\text{O}$ + 2.85 g $\text{Na}_2\text{HPO}_4 \times 2 \text{H}_2\text{O}$ + 9 g NaCl, fill up to 1 l with distilled water/methanol (90/10)
- 0.13 M sodium hydrogen carbonate buffer, pH 8.1:
10.92 g NaHCO_3 fill up to 1 l with distilled water, the pH will be approx. at 8.1

5. Warnings and precautions for the users

Ochratoxin A is a toxic and carcinogenic substance. Particular care should be taken. Avoid contact with the skin (use gloves).

Decontamination of the glassware and toxin-content solutions is best carried out using a sodium hypochlorite (bleach) solution (10 % (v/v)) overnight (adjust solution with HCl to pH 7).

6. Storage instructions

Store the immunoaffinity columns at 2 - 8 °C (36 - 46 °F) - DO NOT FREEZE.

Return any unused columns to their original foil bag and further store at 2 - 8 °C (35 - 46 °F).

No quality guarantee is accepted after the expiration date on the kit label.

7. Preparation of samples

The samples should be stored in a cool place, protected against light.

Ochratoxin A is light sensitive, therefore, avoid exposure of the samples and sample extracts to direct light.

With the following extraction protocol the actual volume of the sample passed through the column is equivalent to 1 g of starting commodity. The detection limit depending on the analytical method can be increased or decreased in the variation of this amount.

7.1. Cereals and malt

- weigh 10 g of ground sample into a screwable vial and add 20 ml of acetonitrile/water (60/40 (v/v))
- mix well for 20 min (horizontal position preferably)
- filter the solution through filter paper by using a frit and keep the filtrate
- for degreasing add 5 ml of n-heptane to 5 ml of the filtrate
- mix for 5 - 10 min
- centrifuge: 5 min / 3000 g / room temperature (20 - 25 °C / 68 - 77 °F)
- remove the upper heptane layer completely
- dilute 2 ml of the degreased filtrate (equivalent to 1 g of sample) with 13 ml PBS buffer
- pass the entire amount (= 15 ml) of the sample solution through the column

Remark:

Further application notes for RIDASCREEN® Ochratoxin A 30/15 (R1311) in combination with RIDA® Ochratoxin A column are available on request. Please contact your local distributor.

- for green and roasted coffee
- for wine
- for dried fruits

Further application notes for RIDASCREEN®FAST Ochratoxin A (R5402) in combination with RIDA® Ochratoxin A column are available on request. Please contact your local distributor.

- for green coffee
- for wine
- for dried fruits

Samples prepared by other methods can also be used for analysis with the RIDA® Ochratoxin A column. However, these samples have to be provided in aqueous solutions containing maximum 10 % methanol.

8. Separation with RIDA® Ochratoxin A column

8.1. Preliminary comments

Bring the columns to room temperature (20 - 25 °C / 68 - 77 °F) before use.

Return unused columns to 2 - 8 °C (35.6 - 46.4 °F), into the original bag.

Remark:

The columns are delivered without an adapter for sample reservoir.

The columns are plugged with caps on top and tip, which have to be removed prior to use. The storage buffer above the gel is rinsed off within the first washing step.

The columns must not dry up during usage. Exceptions are in section 8.2. step 9 (brief removal of residual fluid in the column gel prior to elution) and step 13 (complete removal of residual fluid after elution).

Excessive pressure on the column may cause compression of the gel and consequently low recoveries.

8.2. Clean up

Do not let the columns dry up!

1. equilibrate the column by rinsing with 2 ml PBS/methanol (90/10)
2. fill the column with approx. 1 ml sample extract
3. attach suitable adapter on top of the column and use syringe as sample reservoir
4. fill syringe with the residual sample extract
5. pass sample extract slowly and continuously through the column (flow rate: approx. 1 drop/sec), use positive pressure with syringe or absorb when using vacuum unit
6. discard the passed solution
7. rinse the column with 10 ml PBS/methanol (90/10)
8. discard the passed solution
9. ensure that all liquid is removed from the column by pressing air or N₂ through the column
10. remove syringe and place a clean and closable vial directly below the column
11. elute with 1 ml methanol (100 %); methanol has to pass **slowly** through the column (flow rate: approx. 1 drop/sec) to ensure complete elution of the toxin
12. if the eluent passed too fast (faster than 10 sec), collect the eluate and pass it again through the column
13. press air thoroughly through the column or absorb with vacuum for 30 sec to collect all eluat residues

9. Detection and quantification

9.1. RIDASCREEN® Ochratoxin A 30/15 (R1311)

- evaporate the toxin containing eluate (equivalent to 1 g cereals/malt) to complete dryness at approx. 60 °C (140 °F), if possible by means of a weak nitrogen stream
- redissolve the dried residue in 1 ml of 0.13 M sodium hydrogen carbonate buffer and dissolve thoroughly (vortex or 3 - 5 min in an ultrasonic bath)
- use 50 µl per well in the RIDASCREEN® Ochratoxin A 30/15 assay

Please note:

In order to obtain the ochratoxin A concentration in ng/kg actually contained in a sample the concentration read from the calibration curve has to be multiplied by the corresponding dilution factor.

If 15 ml of cereals/malt extract is passed through the column in accordance with the protocol given in 7. Preparation of samples, the resulting dilution factor is as follows:

cereals/malt: 1

measuring range: between 50 ng/kg (ppt) and 1800 ng/kg (ppt)

remarks for storage of samples and sample extracts:

original samples: always keep refrigerated, dry, protected against light and well sealed

methanolic eluates: store in a cool place protected against light; protected against light and frozen the eluate can be stored for 2 to 3 months

store sample extracts in sodium hydrogen carbonate buffer (see 9.1. after evaporation of the eluate) at 2 - 8 °C (35.6 - 45.4 °F) only for 2 - 3 days

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG