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# Set of 3 processed Gliadin Assay Controls

## 1. Intended Use

The processed assay controls are used to control extraction, test performance and handling of the RIDASCREEN<sup>®</sup> Sandwich ELISA (Art. No. R7001 and R7002) and RIDA<sup>®</sup>QUICK assays (Art. No. R7003 and R7004) for detection of Gliadin. The assay controls need to be extracted with the Cocktail (patented) (R7006/R7016, official R5-Mendez method).

## 2. Content

Each assay control contains 1.5 g homogenized snack sample. The concentration of the respective assay control after extraction with the Cocktail (patented) is specified in the Certificate of Analysis.

For production of the assay controls, snack samples were ground (processed food without tannin or polyphenol containing ingredients), which were well characterized referring their gliadin-content [1]. Furthermore, the defined snack samples were analyzed in an official interlaboratory study of AACCI and PWG (Prolamin Working Group) in 2012. For the preparation of the assay controls, the snack samples were partly remixed.

## 3. Storage

Store assay controls at 2 - 8 °C (35.6 - 46.4 °F) in the dark.

## 4. Preparation of samples

The extraction has to be done as described in the respective test kit instructions. The samples need to be prepared with the Cocktail (patented) (R7006/R7016).

## 5. Test implementation

The further procedure is explained in the respective test kit instructions.

## 6. Results

The characteristic of RIDASCREEN<sup>®</sup> Sandwich ELISA is the quantitative analysis of prolamins from wheat (gliadin), rye (secalin) and barley (hordein) in food declared as gluten-free. A gliadin concentration in mg/kg (ppm) is read out from

the RIDA<sup>®</sup>SOFT Win calibration curve including the respective sample dilution. The quantitative target range for the assays R7001 and R7002 is specified in the certificate provided with the controls.

If the results do not correspond to the target range, a trouble shooting of the assay procedure has to be carried out.

The RIDA<sup>®</sup>QUICK assays provide a qualitative result. According to the test strips' limit of detection for processed samples, the control samples have as target a negative or positive result respectively.

Often airborne cereal dust or dirty laboratory equipment lead to gliadin contamination in the assay. It is recommended to clean surfaces, pipettes, glass vials and mincers before use with 40 % ethanol or propanol.

## Remark

The gliadin assay controls are **not** a reference material. Incurred samples should be used for the validation of the RIDASCREEN<sup>®</sup> Gliadin tests. Applications can be provided by R-Biopharm.

## Literature

[1] Koehler P, Schwalb T, Immer U, *et al.* AACCI Approved Methods Technical Committee Report: Collaborative study on the immunochemical determination of intact gluten using an R5 sandwich ELISA. *Cereal Foods World* 2013; **58**: 36-40.

The data correspond to our present state of technology and provide information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

### R-Biopharm AG

Postanschrift / Postal Address:

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64297 Darmstadt, Germany

Sitz / Corporate Seat: Pfungstadt

Tel.: +49 (0) 61 51 - 81 02-0

Fax: +49 (0) 61 51 - 81 02-40

E-mail: [info@r-biopharm.de](mailto:info@r-biopharm.de)

[www.r-biopharm.com](http://www.r-biopharm.com)

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are registered trademarks of R-Biopharm AG