**SureFood<sup>®</sup> FISH ID**  
**Merluccius merluccius IAAC**  
**(R&D Version)**

Art. No. S6311

50 rxn

User Manual

 **Inhalt**

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 **Content**

1	General Information .....	7
1.1	Description .....	7
1.2	Limit of Detection .....	7
1.3	DNA-preparation .....	7
1.4	Kit components and storage .....	7
1.5	Additionally required equipment and materials .....	7
1.6	Setup.....	7
1.7	Detection channel Set-up .....	8
2	Qualitative Analysis .....	9
2.1	Protocol .....	9
2.1.1	Preparation of the master-mix .....	9
2.1.2	Preparation of the real-time PCR-mix .....	9
2.2	Interpretation of results .....	10
3	Further Information .....	10
3.1	Product Information .....	10
3.2	Technical Support .....	10
3.3	Distribution and Ordering .....	10

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## 1 General Information

### 1.1 Description

The test detects European hake (*Merluccius merluccius*) DNA. Each reaction contains an internal amplification control and an internal detection assay for vertebrates DNA (IAAC). The real-time PCR assay can be used with established real-time PCR instruments, equipped for detection of two fluorescence emissions at 522 nm and 553 nm at the same time.

### 1.2 Limit of Detection

The SureFood® FISH ID Merluccius merluccius IAAC (R&D Version) real-time PCR is developed for the detection of European hake in muscle meat mixture at a relative amount of 1 %. The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

### 1.3 DNA-preparation

For DNA-preparation the use of SureFood® PREP Basic and for highly processed food and feed the use of SureFood® PREP Advanced is recommended.

### 1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	1 x 1100 µl	Yellow
2	Taq Polymerase	1 x 10 µl	Red
3	Positive Control	1 x 100 µl	Light Blue

Store all reagents at –20°C and protected from light.

### 1.5 Additionally required equipment and materials

- DNA-Extraction kit (e.g. SureFood® PREP Basic or SureFood® PREP Advanced)
- real-time PCR instrument with two detection channels (522 nm and 553 nm)
- real-time PCR consumable (plates, tubes, foils, caps)
- pipettes with filter tips
- unpowdered disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

### 1.6 Setup

	<b>Blockcycler / R-Biopharm RIDA®CYCLER</b>	<b>Rotorcycler / LTF MyGo Pro</b>
Initial Denaturation (HOLD)	5 min, 95°C	1 min, 95°C
Cycles	35	35
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

**1.7 Detection channel Set-up**

<b>Real-time PCR device</b>	<b>Detection</b>	<b>Detection channel</b>	<b>Quencher</b>	<b>Note</b>
<b>Agilent Mx3005P</b>	<i>Merluccius merluccius</i>	FAM	+	
	IAAC	HEX	+	
<b>Agilent AriaDx</b>	<i>Merluccius merluccius</i>	FAM	+	
	IAAC	HEX	+	
<b>Applied Biosystems 7500</b>	<i>Merluccius merluccius</i>	FAM	None	Check the passive reference option ROX is none.
	IAAC	VIC	None	
<b>Bio-Rad CFX96</b>	<i>Merluccius merluccius</i>	FAM	+	
	IAAC	VIC	+	
<b>R-Biopharm RIDA®CYCLER</b>	<i>Merluccius merluccius</i>	green	+	
	IAAC	yellow	+	
<b>LTF MyGo Pro</b>	<i>Merluccius merluccius</i>	FAM	+	
	IAAC	VIC	+	
<b>Qiagen Rotor-Gene Q</b>	<i>Merluccius merluccius</i>	green	+	
	IAAC	yellow	+	
<b>Roche LightCycler® 480 II</b>	<i>Merluccius merluccius</i>	465-510	+	
	IAAC	533-580	+	
<b>Roche cobas® z 480 Analyzer</b>	<i>Merluccius merluccius</i>	465-510	+	
	IAAC	540-580	+	

## **2 Qualitative Analysis**

### **2.1 Protocol**

#### **2.1.1 Preparation of the master-mix**

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay as well as for the inhibition control.

Recommended control reactions for the specific PCR assay: negative control, extraction control, positive control. The reaction mix contains an internal amplification control (IAC) per reaction.

Reactions needed for the qualitative *Merluccius merluccius* detection:

3 reactions for controls (1x no-template control, 1x extraction control, 1 positive control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

<b>Components of the master-mix</b>	<b>Amount per reaction</b>	<b>10 reactions (with 10% excess)</b>
Reaction Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

Mix each master-mix well and centrifuge shortly before use.

#### **2.1.2 Preparation of the real-time PCR-mix**

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

## 2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions need to give the correct results.

*Merluccius merluccius* DNA is detected in the FAM-channel. In the VIC-channel it is possible to detect animal DNA in the sample as well as the amplification control in a sample with no animal DNA inside.

A sample is stated **positive** for *Merluccius merluccius*, if the sample DNA shows amplification in the FAM-channel. A sample is stated **negative** for *Merluccius merluccius*, if the sample DNA shows no amplification in the FAM-channel and if the internal control (VIC-channel) of the sample is positive.

If the internal control (VIC-channel) of the sample DNA is detected significant, before the signal of the negative control (master-mix without DNA) the sample contains animal DNA. Is the CP-value of the internal control (VIC-channel) in the range of the negative control the sample contains no PCR-inhibiting substances but only a low amount or no animal DNA.

If the sample DNA and the internal signal are **negative** the sample contains PCR-inhibiting substances. Under these circumstances an evaluation of the samples is not possible. DNA isolation and purification for the sample need to be improved.

**Note:** If the sample contains more than one animal species the DNA mixture can have a competitive influence on the absolute fluorescence. The lower the concentration of the determinant DNA is in a mixture of animal DNAs the lower is the fluorescence level of the amplification curve.

## 3 Further Information

### 3.1 Product Information

- Detailed information about setup of several real-time PCR devices  
(Download: [www.congen.de/en/company/downloads](http://www.congen.de/en/company/downloads))
- Validation Report

### 3.2 Technical Support

For further questions please contact your distributor or send an e-mail to [sales@r-biopharm.de](mailto:sales@r-biopharm.de).

### 3.3 Distribution and Ordering

R-Biopharm AG  
An der neuen Bergstrasse 17,  
64297 Darmstadt, Germany  
Phone: +49 (0) 61 51 - 81 02-0  
Fax: +49 (0) 61 51 - 81 02-20  
E-Mail: [orders@r-biopharm.de](mailto:orders@r-biopharm.de)  
[www.r-biopharm.com](http://www.r-biopharm.com)