RIDASCREEN[®] FAST Soya

Art. No. R7102

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RIDASCREEN[®]FAST Soya

Brief information

RIDASCREEN[®]FAST Soya (Art. No. R7102) is a sandwich enzyme immunoassay for the quantitative analysis of soya proteins in untreated and processed food and beverages.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 48 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation: homogenization and extraction

Time requirement: sample preparation (for 10 samples)approx. 25 min
test implementation (incubation time)30 min

Kalibrator: soya protein (soya flour with 39 % soya protein)

Limit of detection: 0.24 mg/kg (ppm) soya protein

Limit of quantification: 2.5 mg/kg (ppm) soya protein

Specificity: The antibodies specifically detect heated soya proteins.
(The sample extraction at **100 °C for 10 min** ensures that all soya proteins denature and can be detected in the ELISA).

There is a cross-reactivity to legumes of the tribe Phaseoleae (e.g. beans), the genus *Vicia* (e.g. fava bean), to dried peas / pea flour and to peanut.

Cross reactivities of the used antibodies have been determined for the pure food (e.g. corn flour).

In a composed / processed food (e.g. corn bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website www.r-biopharm.com/products/food-feed-analysis.

Related products

RIDA[®]QUICK Soya (Art.No. R7103)

Set of 3 processed Soya Assay Controls (Art.No. R7132)

SureFood[®] ALLERGEN ID 4plex Soya/Celery/Mustard + IAC (Art.No. S3401)

SureFood[®] ALLERGEN ID Soya (Art.No. S3101)

SureFood[®] ALLERGEN QUANT Soya (Art.No. S3201)

1. Intended use

RIDASCREEN[®]FAST Soya (Art. No. R7102) is a sandwich enzyme immunoassay for the quantitative analysis of native and processed soya proteins in food. Such products are beverages and food as sausages, dressings, bakery products, ice, chocolate, soup, sauce, margarine, etc.

2. General

Soya proteins can be present as an ingredient or as a contamination in raw and processed products. According to the **regulation (EU) No. 1169/2011**, soya must be declared on food labels as it can induce allergic reactions. Similar regulations exist e.g. in the USA, Canada, Australia and New Zealand.

Soybeans (mature seeds, raw) contain approx. 35 % proteins. Therefore, soya is one of the few plants which are considered as a valuable source of proteins. Thus, soya is widely used as substitute for proteins from animal sources (e.g. in tofu, soya milk or soya yoghurt). Besides the high protein content, soya beans are also rich in fat (approx. 20 %) and are used for the productions of oils and fats.

The usage of soya as food and feed has strongly increased in the last decades.

3. Test principle

The wells of the microtiter strips are coated with specific antibodies to soya proteins. By adding standards and samples to the wells, any soya protein present will bind to the specific antibodies. In a washing step components not bound are removed. Then antibody conjugated to peroxidase is added. This antibody conjugate is bound to the Ab-Ag-complex. An antibody-antigen-antibody (sandwich) complex is formed. Any unbound conjugate is then removed in a washing step. The detection of soya protein takes place by adding Substrate/Chromogen. The enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is proportional to the soya proteins concentration of the sample. The result is expressed in mg/kg soya proteins.

4. Reagents provided

Each kit contains sufficient materials for 48 measurements (including standard analyses).

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		48 wells
4 x Extractor 3	Brown	Ready to use		30 ml
Allergen extraction buffer	Green	Concentrate	10x	100 ml
Standard 1*	Transparent	Ready to use	0 mg/kg	1.3 ml
Standard 2*	Transparent	Ready to use	2.5 mg/kg	1.3 ml
Standard 3*	Transparent	Ready to use	5.0 mg/kg	1.3 ml
Standard 4*	Transparent	Ready to use	10.0 mg/kg	1.3 ml
Standard 5*	Transparent	Ready to use	20.0 mg/kg	1.3 ml
Wash buffer	Brown	Concentrate	10x	100 ml
Conjugate	Red	Concentrate	11x	0.7 ml
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use		10 ml
Stop solution	Yellow	Ready to use		14 ml

*) The **dilution factor of 100** for the sample has already been considered when labeling. Therefore, the soya protein concentration of samples can directly be read from the standard curve.

5. Reagents required but not provided

5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- centrifuge, centrifugal vials
- shaker
- water bath (60 °C / 140 °F and 100 °C / 212 °F)
- laboratory mincer / grinder, pestle and mortar, Ultra-Turrax or homogenisator
- paper filter
- graduated pipette
- variable 20 - 200 µl and 200 - 1000 µl micropipettes
- if possible multichannel pipette and/or multistep pipette

5.2. Reagents:

- distilled or deionized water
- for the sample preparation of tannin and polyphenol containing food: casein (e.g. Sigma C5890) and polyvinylpyrrolidon (PVP) K15 or K17 (e.g. AppliChem A2258)

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

Reddish substrate/chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- any bluish coloration of the red stained substrate/chromogen prior to test implementation
- a value of less than 0.8 absorbance units ($A_{450\text{ nm}} < 0.8$) for standard 5
- odorless Extractor 3 after repeated opening of the bottle; Extractor 3 should smell like mercaptoethanol (rotten eggs). (For checking, carefully fan the vapor of the open bottle).

9. Preparation of Samples

Working devices such as mill, glass vials or spatulas must be cleaned before and after each sample preparation to remove any remains of allergen and to avoid contamination.

The **Allergen Extraction Buffer** is provided as a **10fold concentrate**. Before dilution, dissolve possibly formed crystals by heating (water bath 37 °C / 98.6 °F) and mix well. Dilute the heated concentrate 1:10 (1+9) with distilled water (e.g. 100 ml concentrate + 900 ml dist. Water). The diluted Allergen Extraction Buffer can be used for approx. 12 weeks at 2 - 8 °C (35.6 – 46.4 °F). It is recommended to work under a chemical hood, because of β -mercaptoethanol in the Extractor 3. The pH of the samples should be set to a neutral value.

9.1. Extraction solid samples (**without** Casein)

- homogenize a representative amount of sample (5 – 50 g)
- weigh 1 g of the homogenized sample, add 2.5 ml Extractor 3 and 17.5 ml preheated (60 °C / 140 °F) and diluted Allergen Extraction Buffer, close the vial, mix vigorously and put the solution for **10 min at 100 °C / 212 °F** in a boiling water bath (the entire sample tube should be covered with boiling water)
- let the sample cool down shortly (e.g. in an ice bath) and centrifuge 10 min at least 2500 g, if possible at 4 °C (39 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged at high speed for 10 min in reaction caps by using a microcentrifuge)
- dilute the extract 1:5 (1+4) with diluted Allergen Extraction Buffer (e.g. 100 μ l sample extract + 400 μ l diluted Allergen Extraction Buffer)
- use 100 μ l extract per well in the assay

9.2. Extraction liquid samples (**without** Casein)

- take 1 ml of the homogenized sample, add 2.5 ml Extractor 3 and 16.5 ml preheated (60 °C / 140 °F) and diluted Allergen Extraction Buffer, close the vial,

- mix vigorously and put the solution for **10 min at 100 °C / 212 °F** in a boiling water bath (the entire sample tube should be covered with boiling water)
- let the sample cool down shortly (e.g. in an ice bath) and centrifuge 10 min at least 2500 g if possible at 4 °C (39 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged at high speed for 10 min in reaction caps by using a microcentrifuge)
- dilute the extract 1:5 (1+4) with diluted Allergen Extraction Buffer (e.g. 100 µl sample extract + 400 µl diluted Allergen Extraction Buffer)
- use 100 µl extract per well in the assay

9.3. Extraction of tannin and polyphenol containing food like spices, chocolate and coffee (**with** Casein and PVP)

To the finally diluted Allergen Extraction Buffer (AEB), add 0.25 % casein and 1 % PVP, this aids extracting soy proteins from different spices, chocolate and other polyphenol containing matrices (see validation data, available at R-Biopharm). Components are better soluble if the solution is preliminarily heated to 60°C / 140 °F.

- homogenize a representative amount of sample (5 – 50 g)
- weigh 1 g of the homogenized sample, add 2.5 ml Extractor 3 and 17.5 ml preheated (60 °C / 140 °F) and diluted AEB (containing 0.25 % casein and 1 % PVP), close the vial, mix vigorously and put it for **10 min at 100 °C / 212 °F** in a boiling water bath (the entire sample tube should be covered with boiling water)
- let the sample cool down shortly (e.g. in an ice bath) and centrifuge 10 min at least 2500 g if possible at 4 °C (39 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged at high speed for 10 min in reaction caps by using a microcentrifuge)
- dilute the extract 1:5 (1+4) with diluted Allergen Extraction Buffer (e.g. 100 µl sample extract + 400 µl diluted Allergen Extraction Buffer)
- use 100 µl extract per well in the assay

The diluted Casein/PVP-Extraction Buffer can be stored at 2 - 8 °C (35 - 46 °F) for 7 days. The sample extracts can be stored at 2 - 8 °C (35 - 46 °F) for 3 days. The extracts can be stored at -20 °C (-4 °F) for several months.

10. Test implementation

10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **conjugate** (bottle with red cap) is provided as a 11fold concentrate. Since the diluted conjugate solution has a limited stability, only the amount which actually is needed should be diluted. Before pipetting, the conjugate concentrate should be shaken carefully. For reconstitution, the conjugate concentrate is diluted 1:11 (1+10) with distilled water (e. g. 200 µl concentrate + 2 ml distilled water, sufficient for 2 microtiter strips).

The **wash buffer** is provided as a 10fold concentrate. Before use, the buffer has to be diluted 1:10 (1+9) with distilled water (i.e. 100 ml buffer concentrate + 900 ml dist. water). Prior to dilution, dissolve eventually formed crystals by incubating the buffer in a water bath at 37 °C (99 °F). The diluted buffer is stable at 2 - 8 °C (35 - 46 °F) for four weeks.

10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than three strips (24 wells) at the same time. In the case of more than three strips, a second uncoated plate (e.g. low binding from Greiner bio-one Cat.-No. 655101) should be used as a pre-plate to avoid a time shift over the microtiter plate. All standards and samples are pipetted into the uncoated plate (at least 150 µl per well) and then quickly transferred to the coated microtiter plate with an 8-channel pipette.

1. Insert a sufficient number of wells into the microwell holder for all standard and samples to be run in duplicate. Record standard and sample positions.
2. Add 100 µl of each standard or prepared sample to separate duplicate wells and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl wash buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
4. Add 100 µl of the diluted conjugate to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).

5. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl wash buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
6. Add 100 µl substrate/chromogen to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 – 77 °F) in the dark.
7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes after addition of stop solution.

11. Results

A special software, the RIDA[®]SOFT Win / RIDA[®]SOFT Win.net (Art. Nr. Z9996), is available for evaluation of the RIDASCREEN[®] enzyme immunoassays. The calculation should be done by use of a cubic spline function. The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit. For quality assurance Soya Assay Controls (R7132) should be used.

In comparison with the certificate, higher values of the absorbance ($A_{450\text{ nm}}$) for the standard curve, especially for the zero standard, may be a result of insufficient washing or allergen contamination.

A further dilution and new detection of the samples is recommended for absorbance values ($A_{450\text{ nm}}$) > standard 5.

Please note

When working in accordance with the regulation stated, the dilution factor is 100. The allergen concentration can be read directly from the standard curve - the sample dilution factor of 100 is already taken into account for the standard concentrations (see 4. *).

For sample dilutions of more than 1:100 the further dilution factor must be considered for the calculation of the soya protein concentration.

E.g. if a sample is tested for 1 % soya flour concentration (corresponding to 10000 mg/kg soya flour or approx. 3900 mg/kg soya protein), the sample has to be diluted 1:70000 (1:100 according to the sample preparation and additional 1:700 with Allergen Extraction Buffer).

In general

Samples tested negative still could contain an allergen contamination below the limit of detection of the assay, or they might contain other allergen components like lipids for example.

Due to the multitude of food types, matrix effects cannot be excluded. In processed food (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented, this may have an impact on the recovery/cross reactivity.

For evaluation of the cross reactivity only one exemplary sample was analyzed, other samples may show a different result, all cross reactivities and exemplary analyzed matrices are described in the validation report.

Recommendation

In order to ensure a high analytical performance we would recommend:

- to analyze each sample material should in duplicates
- to use also allergen-free and allergen- containing (spiked) samples as test controls
- due to the multitude of food types, matrix effects cannot be excluded. To ensure an accurate result spike experiments are recommended
- to perform SureFood[®] PCR for confirmation of the result
- For details using the ChemWell[®] or GEMINI automation please contact sales@r-biopharm.de.

For further product information and application notes, please contact sales@r-biopharm.de.

The product information folder with further information is available on request from your local distributor or R-Biopharm AG.

Further application notes

- ELISA - Swabbing Allergens
- Sample preparation caseinates

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG