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Order No.: QPP1HR 0050

(Version 11/16)

1. Examples of successful detected microorganisms

- *Lactobacillus* spp.
 - L. acetotolerans*
 - L. acidophilus*
 - L. amylolyticus*
 - L. amylovorus*
 - L. backii*
 - L. brevis*
 - L. brevisimilis*
 - L. buchneri*
 - L. casei* / *paracasei*
 - L. collinoides* / *paracasei*
 - L. coryniformis* ssp. *coryniformis*
 - L. coryniformis* ssp. *torquens*
 - L. curvatus*
 - L. delbrückii* ssp. *delbrückii*/ *lactis*
 - L. harbinensis*
 - L. helveticus*
 - L. johnsonii*
 - L. lindneri*
 - L. malefermentans*
 - L. parabrevis*
 - L. parabuchneri* (*frigidus*)
 - L. paucivorans*
 - L. pentosus*
 - L. perolens*
 - L. plantarum*
 - L. paraplantarum*
 - L. rossiae*
- *Lactococcus lactis*
- *Megasphaera* spp.
- *Micrococcus kristinae*
- *Pectinatus* spp.
 - Pectinatus cerevisiiphilus*
 - Pectinatus frisingensis*
 - Pectinatus haikarae*
 - Pectinatus* ssp.
- *Pediococcus* spp.
 - Pediococcus acidilactici*
 - Pediococcus damnosus*
 - Pediococcus dextrinicus*
 - Pediococcus inoptinatus*
 - Pediococcus pentosaceus*
 - Pediococcus parvulus*
 - Pediococcus claussenii*

2. Intended use

Qualitative detection of bacterial contaminations and the hop resistance genes *hor A* and *hor C* as genetic markers for beer-spoilage microorganisms in beer and beer mixes.

3. Test principle

The TaqMan[®] real-time PCR is based on hot-start-PCR and sequence-specific dual labelled probes (FAM/DQ; HEX- or JOE/DQ; 610 or ROX/DQ) which, when accurately hybridised, emit a measurable fluorescent signal of a defined wavelength in the extension phase. The increase of signal is continuously measured in a real-time PCR detection instrument. To avoid false negative PCR-results an Inhibition Control (660 or CY5/DQ) is amplified together in one reaction vessel with the specific sequence. The system contains dUTP. Optional: the use of Uracil-N-Glycosylase will eliminate any contamination with Uracil containing amplicons from former PCRs (the enzyme is not part of this kit).

4. Kit contents

The kit contains sufficient reagents for 50 reactions:

1 x Premix QPP1HR	white cap
1 x QPP1HR mix (freeze-dried, incl. IC-DNA)	dark vial, red cap
1 x ddH ₂ O	colourless cap
1 x Control-DNA (freeze-dried)	yellow cap

5. Storage conditions

The QPP1HR mix and the Control-DNA are freeze-dried, they have to be solved in ddH₂O prior to use (see 6.1).

Do **not** freeze the lyophilized PCR mix and lyophilized Control-DNA.

All PCR reagents, except the Premix should be stored at 4 – 8 °C (39 – 46 °F).

Keep Premix for storage at - 20 °C (- 4 °F). Avoid loss of sensitivity by repeating freezing and thawing more than 3 times. For irregular use aliquot the Premix.

The QPP1HR mix contains the fluorescent labelled probes and should be handled light protected.

All reagents are stable for 6 months, if they are stored correctly

6. Materials required, but not provided

6.1. Instruments

LightCycler™ 480 or real-time PCR machine with channels FAM, HEX (JOE), ROX and CY5

Standard benchtop mini-centrifuge (1.5 – 2.0 mL vials)

Centrifuge for multiwell plates or stripes

Pipettes

“Vortex”

6.2. Reagents and plastic ware

sterile ddH₂O

sterile reaction vessels 1.5 – 2.0 mL

sterile plates or optical tubes (plastic)

sterile filter tips

optional: Uracil N-Glycosylase (0.01 U/μL added to the PCR reaction mix)

Colour Compensation Set for LC480 (GEN-IAL®, PP1T CC LC480 0005)

Use a Colour compensation prior to first use of QPP1HR Kit!

7. PCR

7.1. PCR-Setup

Preparations before first use:

When using the kit for the first time, the freeze-dried kit components have to be shortly centrifuged and carefully resolved:

- add 80 μL sterile ddH₂O to freeze-dried QPP1HR mix
- add 55 μL sterile ddH₂O to the freeze-dried Control-DNA
- after 15 minutes mix well

Before every use thoroughly mix all PCR components and centrifuge briefly.

PCR-reaction Setup:

PCR-Components	volume (μL)
Premix QPP1HR	13.5
QPP1HR mix	1.5
Sample-DNA	5.0*
Total volume	20.0

* if using the Simplex easy DNA kit: add 2.5 μL DNA

1. Prepare a mastermix by mixing Premix and QPP1HR mix.
2. Multiply said volumes with the number of PCR preparations including controls (positive control, negative control, extraction control), taking into account pipette reserves of approximately 5-10 %.
3. Divide 15 μL of the PCR-mastermix among the individual reaction vessels, making sure that, prior to the first filling, the tip of the pipette has been moistened.
4. Add the sample DNA using a fresh tip with each DNA filling and mix it with the mastermix at the same time, (with DNA volume less than 5 μL add ddH₂O to the end volume of 20 μL per reaction). Add 5 μL of the Control-DNA, 5 μL extraction control and 5 μL of ddH₂O for the negative control reaction.
5. Close the tubes immediately and centrifuge them according to the recommendations of the real-time PCR instrument manual.
6. Place tubes/plate/capillaries into the thermo block/carousel of the instrument, close lid and start run.

Work swiftly to avoid warming up and keep away from light.

7.2 PCR-Program (LightCycler 480 Software release 1.5.0 SP1):

7.2.1 PCR-Program LC480

1. Click the button **tool** on the right side in the window **LightCycler 480 Software release 1.5.0. SP1**
2. Click the button **Detection formats** at the left side of the menu bar
3. Click **New** in the window **Detection formats** and name the experiment, for example: QPP1HR
4. Open the window **Filter Combination Selection** and choose the following filter combinations: 465-510 / 533-580 / 533-610 / 618-660.
5. Open the window **Selected Filter Combination** list and add the following amounts

Excitation Filter	Emission Filter	Name	Melt factor	Quant fact	Max. Integ. Time
465	510	465-510 Fam	1	10	2
533	580	533-580 Hex	1	10	2
533	610	533-610 610	1	10	2
618	660	618-660 660	1	10	2

6. **Close** the window
7. Click the button **New Experiment** on the right side of the menu bar
8. From the pull-down menu **Detection formats** choose the defined experiment, click the button **Customize** and check the detection formats. All of them have to be activated
9. Click the button **ok**
10. Define the following program:

1. Program Name: **Heat**

Cycles **1** Analysis Mode **None**

Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)	Sec Target (°C)	Step size (°C)	Step Delay (cycles)
95	None	00:15:00	4.40		0	0	0

2. Program Name: **Ampli**

Cycles **35** Analysis Mode **Quantification**

Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)	Sec Target (°C)	Step size (°C)	Step Delay (cycles)
95	None	00:00:10	4.40		0	0	0
60	Single	00:00:20	2.20		0	0	0

3. Program name : **Cool**

Cycles **1** Analysis Mode **None**

Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)	Sec Target (°C)	Step size (°C)	Step Delay (cycles)
40	None	00:00:20	2.20		0	0	0

Optional: Saving the program as **run template**. The template button allows to select and apply a template to the currently open object and to save the currently open object as a template. Click the clamp beside the button **Apply Template** . Click **Save as template** and save the file in **templates**.

11. Click the button **Subset editor** on the left side

12. Click the button **+** and **New Subset 1** appears

13. Mark the wells in the **New Subset 1 Settings** window

14. Click the button **Apply**

15. Click the button **Sample editor** on the left side

16. **Very important:** Activate In the window **Step 1 Select Workflow: Abs.Quant**

17. Choose the subset **New Subset** in the window **Step 2 Select Samples**

18. Define your probes in the **Sample table**

7.2.2 PCR program for other real-time machines

The Premix contains no ROX. This must be considered according to the settings of the real-time machine. **ABI 7500:** Under "Assign Targets and Samples" in "Select the dye to use as the passive reference" choose „none“.

For the use of UNG the thermal cycler program has to be changed according to manufacturers` instructions.

Step	Time	Temp.	
Initial denaturation of DNA	15 min	95 °C	
Cycling Denaturation	15 sec	95 °C	Cycle 35 x
Cycling Annealing/ Elongation	20 sec	60 °C	
Cooling	15 sec	40 °C	

8. Evaluation

The evaluation has to be made according to the data analysis program recommended by the real-time instrument manufacturer.

For LC480: Activate Colour Compensation

***Lactobacillus/ Pediococcus*-DNA:** FAM-channel (LC480: 465-510 nm)

***Megasphaera/ Pectinatus*-DNA:** ROX-channel (LC480: 533-610 nm)

Hor A and Hor C: HEX-channel (LC480: 533-580 nm)

Inhibition Control-DNA: CY5-channel (LC480: 618-660 nm)

A sample is ***Lactobacillus/Pediodoccus* positive**, if there is a detectable fluorescence increase in the **FAM-channel (465-510 nm)** and the negative controls show no amplification. The positive controls should have a positive fluorescence signal. The Inhibition Control in the CY5-channel (618-660 nm) may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). In the negative controls it has to be positive.


A sample is ***Megasphaera/ Pectinatus* positive**, if there is a detectable fluorescence increase in the **ROX-channel (533-610 nm)** and the negative controls show no amplification. The positive controls should have a positive fluorescence signal. The Inhibition Control in the CY5-channel (618-660 nm) may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). In the negative controls it has to be positive.

A sample is **Hor A and/or Hor C positive**, if there is a detectable fluorescence increase in the **HEX-channel (533-580 nm)** and the negative controls show no amplification. The positive controls should have a positive fluorescence signal. The Inhibition Control in the CY5-channel (618-660 nm) may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). In the negative controls has to be positive.

A sample is **negative**, if there is no detectable fluorescence increase in the different channels and the positive controls have a positive fluorescence signal. The negative controls show no amplification. The Inhibition Control in the CY5-channel (618-660 nm) has to be positive in the sample and in the negative controls, a false negative result due to inhibitory effects is then excluded.

analysis flowchart

LC480: analysis after activating Colour Compensation

Note:

The polymerase-chain reaction (PCR) is protected by patents and requires a licence from Hoffmann-LaRoche Inc..

The provided product does not authorise the purchaser for the commercial use of this method.

GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.