**SureFast® Hepatitis E PLUS (100 Reakt.)**

Art. Nr. F7142

Januar 2017

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**Description**

SureFast® Hepatitis E PLUS is a real-time RT-PCR for the direct and qualitative detection of hepatitis E virus. The test contains an Internal Control RNA (ICR, consisting of MS2-bacteriophage) as an internal control of sample preparation procedure and to determine possible PCR-inhibition. The real-time PCR assay can be used with established real-time PCR instruments, equipped for detection of two fluorescence emissions at 522 nm and 553 nm (FAM and VIC/HEX) at the same time. The technical validation of instruments was performed on Agilent Mx3005P, BioRad CFX 96, Roche LightCycler® 480 II, Roche Cobas z 480 Analyzer, Roche LightCycler® 2.0<sup>2</sup>, Applied Biosystems 7500, LTF MyGo Pro and Qiagen RotorGene Q.

**Limit of Detection**

The SureFast® Hepatitis E PLUS real-time RT-PCR has a limit of detection of  $\leq 25$  RNA copies. The assay limit of detection depends on RNA preparation and RNA content.

**RNA-preparation**

For RNA-preparation the use of SureFast® PREP DNA/RNA Virus is recommended.

**Kit components and storage**

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 700 µl	Yellow
2	PP-Mix	1 x 770 µl	Light Green
3	Enzyme Mix	1 x 80 µl	Red
R	Internal Control RNA	2 x 1800 µl	Brown
N	PCR Water	1 x 500 µl	White
P	Positive Control	1 x 100 µl	Light Blue

Store all reagents at -20°C and protected from light.

**Additionally required equipment and materials**

- RNA-Extraction kit (e.g. SureFast® PREP DNA/RNA Virus)
- Real-time PCR instrument, equipped with two detection channels 522 nm and 553 nm
- Real-time PCR consumables (plates, tubes, foils, caps)
- Pipettes with filter tips
- disposal gloves
- Vortexmixer
- Centrifuge with a rotor for the reaction tubes

**Protocol****1. RNA-preparation**

For RNA-preparation the use of SureFast® PREP DNA/RNA Virus is recommended.

The test assay contains an Internal Control RNA (ICR), which can either be used as PCR inhibition control or as extraction control for the sample preparation procedure and as a PCR inhibition control.

If the ICR is used only as a PCR inhibition control, 1 µl of the ICR should be added to the Master-Mix.

If the ICR is used as an extraction control for the sample preparation procedure and as PCR inhibition control, 20 µl of the ICR should be added during extraction procedure. The ICR should always be added to the specimen-lysis buffer mixture and must not be added directly to the specimen.

<sup>2</sup>For the use of the Roche LightCycler® 2.0 a Color Compensation is necessary. The SureCC Color Compensation Kit II (Art. No. F4010) must be used for the color compensation of such devices.

## 2. Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions). Recommended control reactions: negative control, positive control and blank of extraction. The test assay contains an Internal Control RNA (ICR), which can either be used as PCR inhibition control or as extraction control.

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions for ICR as extraction and PCR inhibition control:

Components for master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	12.5 µl	137.5 µl
PP-Mix	6.9 µl	75.9 µl
Enzyme Mix	0.7 µl	7.7 µl
<b>Total volume</b>	<b>20.1 µl</b>	<b>221.1 µl</b>

Mix each master-mix well and centrifuge shortly before use.

Example for the calculation and preparation of 10 reactions for ICR only as PCR inhibition control:

Components for master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	12.5 µl	137.5 µl
PP-Mix	6.9 µl	75.9 µl
Internal Control RNA	1.0 µl	11.0 µl
Enzyme Mix	0.7 µl	7.7 µl
<b>Total volume</b>	<b>21.1 µl</b>	<b>232.1 µl</b>

Mix each master-mix well and centrifuge shortly before use.

## 3. Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes or wells.
- For the negative control pipette 5 µl of PCR Water into the designated tubes and close them.
- Pipette 5 µl of sample RNA into the designated tubes and close them.
- Pipette 5 µl of the Positive Control into the designated tubes and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

4. Setup

	<b>Rotorcyclor/ LightCycler® 480</b>	<b>Rotorgene/ Blockcyclor</b>
Reverse Transkription	10 min, 58 °C	10 min, 58 °C
Initial Denaturation (HOLD)	1 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	10 sec, 95°C	15 sec, 95°C
Annealing/Extension (CYCLE)	15 sec, 60°C	30 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	
Fluorescence Detection Setup	Detection: End of extension phase  Detection system: diverse devices            FAM-channel, Quencher: BHQ Rotor-Gene Q                Green LC480 II                        465 nm, 510 nm  Internal Control RNA: diverse devices            VIC/HEX-channel, Quencher: BHQ Rotor-Gene Q                Yellow LC480 II                        533 nm, 580 nm	
Detailed information on the setup of several real-time PCR devices is available at the CONGEN homepage: <a href="http://www.congen.de/en/company/downloads/">http://www.congen.de/en/company/downloads/</a>		

**Interpretation of results**

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer. The control reactions need to give the correct results.

A sample is stated **positive** for hepatitis E, if the sample RNA shows an amplification in the FAM-channel. A sample is stated **negative** for hepatitis E, if the sample RNA shows no amplification in the FAM-channel and the internal amplification control (VIC/HEX-channel) of the sample is **positive**.

If the sample RNA and the internal amplification control are **negative** the sample contains PCR-inhibiting substances or the sample preparation was not successful. Under these circumstances an evaluation of the samples is not possible. RNA isolation and purification for the sample need to be improved.

**Product Information**

- Validation Report

**Technical Support**

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de).