RIDA® QUICK Aflatoxin RQS ECO

REF R5209

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RIDA®QUICK Aflatoxin RQS ECO

Brief information

RIDA®QUICK Aflatoxin RQS ECO (R5209) is an immunochromatographic test for the quantitative determination of aflatoxin in corn (maize). The test uses an aqueous extraction method.

The test kit contains 20 test strips, each of which can be used for one analysis. All reagents required for the assay are contained in the test kit. Results are evaluated with the RIDA®SMART APP software (ZRSAM1000) or the RIDA®QUICK SCAN reader.

Sample preparation: homogenization and extraction

Time requirement: sample preparation (for 10 samples) approx. 10 min
test implementation (incubation time) 5 min

Limit of detection (LOD): < 2 µg/kg (ppb)

Limit of quantification (LOQ): 2.89 µg/kg

Measuring range: 2 - 75 µg/kg and 50 - 300 µg/kg
Quantitative evaluation with the RIDA®SMART APP

Specificity: The RIDA®QUICK Aflatoxin RQS ECO test detects aflatoxin in corn. Further information is contained in the validation report.

The specificity of the RIDA®QUICK Aflatoxin RQS ECO test was established by determining the cross-reactivity to the corresponding substances in the buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

Related products

RIDASCREEN® Aflatoxin Total (R4701)
RIDASCREEN®FAST Aflatoxin (R5202)
RIDA®QUICK Aflatoxin RQS (R5208)
RIDA®SMART APP Software (ZRSAM1000)
RIDA® Aflatoxin column (R5001 / R5002)
Trilogy® Certified reference material (ISO 17034)

1. Intended use

RIDA®QUICK Aflatoxin RQS ECO is an immunochromatographic test in strip format for the quantitative determination of aflatoxin in corn (maize). The test uses an aqueous extraction method.

2. General information

Aflatoxins are secondary metabolites of the fungi species *Aspergillus flavus*, *parasiticus* and *nomius*. These fungal species are found in humid, tropical areas and can contaminate plant-based food and feed. This contamination usually takes place in the countries of cultivation. Aflatoxins belong to the strongest natural occurring carcinogenic substances.

3. Test principle

The immunochromatographic test in the form of a test strip is based on an antigen-antibody reaction. A specific anti-aflatoxin antibody detects aflatoxin in the sample. During incubation of the test strip, a band pattern (test line / control line) forms that is used to determine the concentration of aflatoxin. The intensity of the test line depends on the aflatoxin concentration of the sample. It increases as the aflatoxin concentration increases. It must be possible to detect the control line after the reaction has taken place in order to verify the functioning of the test. The control line weakens as the aflatoxin concentration in the sample increases.

The test strip is evaluated using the RIDA®SMART APP software (ZRSAM1000) or the RIDA®QUICK SCAN reader.

4. Reagents provided

Each kit contains sufficient materials for max. 20 determinations. Each test kit contains:

| Component | Cap color | Format | Volume |
|--------------------------|-------------|-----------------|---|
| 20 x Test strip | - | ready to use | One determination per test strip (separately packed) |
| 1 x Extraction solution | green | concentrate 10x | 45 ml |
| 1 x Dilution buffer | transparent | ready to use | 20 ml |
| 1 x RIDA®SMART APP cover | - | ready to use | Lot-specific evaluation cover for use with RIDA®SMART APP |

5. Equipment and Reagents required but not provided

5.1 Equipment

- Grinder
- Weighing scale
- Graduated cylinder
- 1 L bottle
- Shaker (e.g. Rock-it 360, Trilogy® Analytical Laboratory Inc. or equivalent)
- Centrifuge (≥ 2000 g)
- 100 μ l pipette (e.g. R-Biopharm FP 100, Art. No. Z0007)
- Latest version of the RIDA®SMART APP software (ZRSAM1000) activated and installed on a compatible smartphone. Alternatively RIDA®QUICK SCAN (Art. No. ZG5005) with software version V1.0.1.7 2.0 or higher.

5.2 Reagents

- Distilled or deionized water

6. Warnings and precautions for the users

This test should be carried out only by trained laboratory employees. The instruction for use must be strictly followed.

The test strips are sensitive to moisture. Moist test strips can have a negative impact on the test result. For this reason, it is important to protect the strips against moisture. Do not remove the test strips from the test strip packaging until immediately prior to use in the test.

This kit may contain hazardous substances. Please refer to the component safety information in the material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

Ensure the proper and responsible disposal of all reagents and materials after their use. For disposal, please adhere to national regulations.

7. Storage instructions

Store the kit at 2 - 8 °C (36 - 47 °F). Do not freeze the test strips or any test kit components.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Sample preparation

To obtain the ready-to-use extraction buffer, transfer 1 part of the extraction solution (concentrate 10x) into 9 parts of distilled or deionized water and homogenize for 5 min. The ready-to-use extraction buffer can be stored for approx. 8 - 10 weeks at 2 - 8 °C. If not all test strips are used within 8 - 10 weeks it is recommended to prepare small volumes of the ready-to-use extraction buffer.

Bring all components necessary for sample preparation to room temperature (20 - 25 °C / 68 - 77 °F) before use in the test and perform sample preparation at room temperature.

The sample should be stored in a cool place, protected against light.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

8.1 Measurement range 2 - 75 µg/kg

- Weigh 10 g of ground and homogenized sample and add 20 ml of ready-to-use extraction buffer
- Shake for 5 min (on shaker)
- Centrifuge (3 min at 2000 g) until a particle-free solution can be obtained
- Use 100 µl of the particle-free solution in the test

8.2 Measurement range 50 - 300 µg/kg

For an extended measurement range of 50 - 300 µg/kg, further dilute the particle-free solution obtained from sample preparation 8.1.

- Dilute 100 µl of the particle-free solution with 700 µl dilution buffer and mix carefully (= Dilution I)
- Use 100 µl of Dilution I in the test

9. Test procedure

Bring all components necessary for sample preparation to room temperature (20 - 25 °C / 68 - 77 °F) before use in the test and perform sample preparation at room temperature.

- Remove a test strip from its foil pouch
- Apply 100 µl of prepared sample solution (particle-free solution or Dilution I, see 8. Sample preparation) to the application area of the test strip (see Fig. 1)
- Incubate the test strip at room temperature
- After exactly 5 min of incubation, evaluate the test strip with the RIDA®SMART APP or the RIDA®QUICK SCAN

10. Evaluation

The left line in the reaction area is the test band (test line, see Fig. 1). Its appearance and intensity depends on the aflatoxin concentration of the sample. The sample is contaminated with aflatoxin, if the control band (control line) is visible and the test band (test line) is visible.

The right line in the reaction area is a control band (control line) and must be present after each test procedure. If this line is not visible, the test was not performed correctly or the reagents were not in order. In this case, the test should be repeated with a new test strip. If the right band is missing again, please inform your local distributor.

Please note: If the sample was prepared according to 8.1 and is highly contaminated with aflatoxin, the control line will be only faintly visible. In this case the test result might be reported as „invalid result“ (RIDA®SMART APP) or „invalid“ (RIDA®QUICK SCAN reader). In this instance, dilute the sample as described in 8.2 and repeat the test with a new test strip.

To obtain comparable results, the test strips must always be evaluated after 5 min.


Fig. 1: RIDA®QUICK Aflatoxin RQS test strip.

10.1 Evaluation with RIDA®SMART APP (ZRSAM1000)

First, please read the instructions for use of the RIDA®SMART APP software attentively. The test strip evaluation has to be performed by using the RIDA®SMART APP software according to these guidelines. Please find a short description (Quick Starter Guide) of how to use the RIDA®SMART APP on our website www.r-biopharm.com.

For the test strip evaluation the RIDA®SMART APP cover is needed which is included in the test kit. The RIDA®SMART APP cover contains lot-specific information and must be used only with the appropriate lot of the test kit. For further information please check the Quality Assurance Certificate enclosed in the kit.

Please make sure the correct RIDA®SMART APP evaluation cover is present prior to starting the test procedure and test strip evaluation with the RIDA®SMART APP software.

RIDA®SMART APP software applications for evaluation:

| Matrix | Measuring range | RIDA®SMART APP application |
|---------------|------------------------|-----------------------------------|
| Corn (maize) | 2 - 75 µg/kg | 2 - 75 ppb |
| | 50 – 300 µg/kg | 50 – 300 ppb |

10.2 Evaluation with RIDA®QUICK SCAN

First, carefully read the user guide for the RIDA®QUICK SCAN and proceed as instructed. Point 3 contains the description for evaluating the test strips.

To obtain optimal results every test kit contains lot-specific standard curve parameters. These parameters can be found on the Quality Assurance Certificate and have to be used only with this specific lot. The Quality Assurance Certificate contains parameters for several standard curves. Depending on the extraction method used, select the corresponding parameters.

To enter the parameters manually into the RIDA®QUICK SCAN reader use the plain text information on the Quality Assurance Certificate. To enter the data with the barcode scanner use the 2D barcodes. In both cases, the checksum shown on the RIDA®QUICK SCAN screen needs to be identical to the respective checksum on the certificate (number below the 2D barcode).

RIDA®QUICK SCAN evaluation methods:

| Measuring range | RIDA®QUICK SCAN method |
|------------------------|-------------------------------|
| 2 - 75 µg/kg | Afla ECO Low |
| 50 - 300 µg/kg | Afla ECO High |

11. Sensitivity

This test determines aflatoxin contaminations in the measuring range of 2 - 300 µg/kg in combination with the RIDA®SMART APP software or the RIDA®QUICK SCAN reader.

Further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

| Version number | Chapter and title |
|----------------|---|
| 2019-11-25 | Release version |
| 2020-12-10 | Brief information and Related products, 4. Reagents provided, 6. Warnings and precautions for the users, 7. Storage instructions, 9. Test procedure, 10. Evaluation, 10.2 Evaluation with RIDA®QUICK SCAN |

Explanation of symbols

- General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM)

Manufacturer + address

Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.