Premi®Test

RBP 31/02-04/11
ALTERNATIVE ANALYTICAL
METHODS FOR AGRIBUSINESS
Certified by AFNOR Certification
End of validity: 30th August 2022
http://nf-validation.afnor.org

Microbial Screening test for the presumptive detection of antibiotic and sulfonamide residues in particularly fresh meat. NF VALIDATION certification scope includes beef, pork and poultry meat (with the exception of ground meat). The following selected antibiotics were investigated: Penicillins: Penicillin G (CCβ = 6 µg/kg); Amoxicillin (CCβ = 11 µg/kg); Cloxacillin (CCβ = 150 µg/kg)
Tetracyclines : Chlorotetracycline (CCβ = 160 µg/kg); Oxytetracycline (CCβ = 160 µg/kg); Doxycycline (CCβ = 100 µg/kg)
Sulfonamides: Sulfadimethoxine (CCβ = 75 µg/kg); Sulfadiazine (CCβ = 90 µg/kg)
Macrolides: Erythromycin A (CCβ = 200 µg/kg); Tylosine A (CCβ = 90 µg/kg)
The CCβ values were determined by spiking in meat juice; the detection capabilities in other samples may differ.
Cefalosporins, Aminoglycosides, Amphenicols, Quinolones, Tiamulin, and Lincomycin are detected above the Maximum Residue Level (MRL).

2019/03/21 - insert 1011

| |
|--|
| |
|--|

Contents

25 ampoules with *Bacillus stearothermophilus var. calidolactis* in a solid agar medium, disposable pipette tips, perforated cover foil + syringe.

Product description

Premi®Test is a broad spectrum microbial screening test especially developed for the detection of antimicrobial substances, such as antibiotic and sulfonamide residues particularly in fresh meat, in the range of MRL.

Test principle

Premi®Test is based on the inhibition of the growth of *Bacillus stearothermophilus*, a micro-organism very sensitive to many antibiotic and sulfonamide residues.

A standardised number of spores are imbedded in an agar medium with selected nutrients. When meat fluid is added to the Premi®Test and heated at 64 °C, the spores will germinate. The germinated spores will multiply and form an acid when no inhibitory substances are present. This will be visible by a colour change from purple to yellow of the indicator in the ampoule. When antimicrobial residues are present in sufficient amount (above the detection level) no growth will occur and the colour will remain purple.

Caution

This test is sensitive to certain antibiotics and other inhibitory substances, any contamination with such materials should be avoided at all time. It is advised to wash hands thoroughly before starting the test procedure. Use tissue paper or a clean towel to dry hands.

Instruction for use

- Wash hands thoroughly before starting the test procedure.
- Cut out the required number of ampoules without damaging the aluminum foil of the adjacent ampoules.
- Remove foil carefully from the required amount of ampoules (do not open more ampoules as needed).
- Take approximately ± 2 mg of lean meat and use a meat press to extract approximately 250 µl of meat fluid. Alternatively, the meat juice can be obtained with the help of the Multipress or by freeze/thawing the meat. (A technical bulletin on the Multipress is available at www.r-biopharm.com).
- Use a new disposable tip on the syringe for each sample.
- Pipette 100 µl of the fluid onto the agar in the ampoule. Do not distort the agar.
- Allow to stand at room temperature for 20 minutes for a pre-incubation.
- Flush the meat juice away by gently filling and emptying the test ampoule twice with demineralized (demi) water. Wash the ampoules with demineralized water only and do not use tap-water!
- Remove the last water carefully from the test ampoule.
- Close the test ampoule with the foil supplied to avoid evaporation.
- Incubate the test ampoule in a Premi®Test incubator or water bath (64 °C ± 1 °C).
- It is required to use, in parallel, a negative control of the same matrix that is being tested. Read the test results after the negative control has changed color.
- All necessary equipment (meat press, scissors, incubator, and timer) are available in the Premi®Test Starter Kit.
- Um to facilitate the use of the Premi®Test there is a Powerpoint document available on request for all users.

Reading test results

- Read the color only from the bottom 2/3 part of the ampoule.
- A clear color change (from purple towards yellow) indicates the absence of antibiotics / sulfonamides above the detection limit.
- No clear color change indicates the presence of antibiotics and/or sulfonamides at/ or above the limit of detection of the test.

Negative control

The use of a negative control is required (obligatory for NF VALIDATION). Start to check the color of the negative control after an incubation time of 2 h 40 min and then in intervals of 5 min until the color of the negative control has changed from purple to yellow. At that time, the samples should be compared to the negative control, because a delayed measuring (5 min.) can lead to increased values of the detection capability (CCβ). If there is no color change of the negative control after 4 h, repeat the test. The negative control can be prepared by freezing and storing meat juice from a sample tested negative. Never use water as a negative control!

Positive control

It is strongly recommended to test a positive control (protocol is available from R-Biopharm AG) to verify the correct application of the test.

Storage

The ampoules should be stored cool (3 – 10 °C). Attention: DO NOT FREEZE!

Limited Liability

Premi®Test is a screening test and as such 100% accuracy of the test results cannot be guaranteed. Besides, the assessment of colour, in particular that of a yellow/purple result, may differ from person to person. In cases where severe consequences are involved for the user, test results should be confirmed by a validated comprehensive analytical method. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality.

Recommendations for sample-preparations of fish, shrimps, eggs, kidney, liver, urine, blood and feed are available at R-Biopharm AG. These matrices require individual validation by the customer. Supportive information are available at R-Biopharm upon request.

| | | |
|--|---|---|
| R-Biopharm AG Postanschrift / Postal Address: An der neuen Bergstraße 17 64297 Darmstadt, Germany Sitz / Corporate Seat: Pfungstadt Tel.: +49 (0) 61 51 - 81 02-0 Fax: +49 (0) 61 51 - 81 02-40 E-mail: info@r-biopharm.de www.r-biopharm.com | Aufsichtsrat/Supervisory Board: Dr. Ralf M. Dreher (Vorsitzender/Chairman) Vorstand/Board of Management: Christian Dreher (Vorsitzender/Chairman), Dr. Hans FrickeI, Jochen Hirsch, Dr. Peter Schubert Handelsregister/ Commercial Register: Amtsgericht Darmstadt, HRB 8321 Sitz der Gesellschaft/ Corporate Seat: Pfungstadt USt.ID-Nr. / VAT-No.: DE 111 657 409 |  |
|--|---|---|

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| | |
|--|--|
| | |
| | |

1. The Intended Use

Premi[®]Test is a Microbial inhibition test validated for the detection of residues of PENICILLIN in fresh bovine meat.

2. The Intended User

The test kit is designed for use by quality control personnel and others familiar with the antibiotic screening testing. All users of the Premi[®]Test should be trained by a R-Biopharm AG representative or by somebody, who already successfully completed the training. To facilitate the use of the Premi[®]Test there is a PowerPoint document available for all users.

3. Contents

- 25 ampoules containing an agar medium, pH indicator and spores of *Bacillus Stearothermophilus*
- 0.1 ml syringe + 25 disposable tips
- Premi[®]Test covering foil

4. Materials required but not included

- Premi[®]Test Starter Kit, containing the required equipment to do Premi[®]Test:
- R-Biopharm AG heating block (Premi[®]Test Incubator), meat press, pair of scissors, thermometer and timer.
- Instead of the Premi[®]Test Incubator a water bath at 64 °C ± 1.0 °C can be used.
- In case a water bath is used: a rack or floater is required to accommodate the ampoules during incubation.
- Calibrated thermometer appropriate for designated use.
- Refrigerator at 3 – 10 °C (37 - 50 °F) for storage of the Premi[®]Test-ampoules.
- Negative control: muscle tissue juice from a known antibiotic-free sample, previously tested as negative.
- Premi[®]Scan, software to interpret the results of the Premi[®]Test, an HP Scanjet 7400 or 8200 scanner, templates, black cloth and Kodak Q-60 Color chart..

5. Test Performance

- 5.1 Detection Sensitivity (dose response information Premi[®]Test):** 10 ppb for Penicillin G The Premi[®]Test has been validated to detect penicillin G residues in fresh bovine muscle tissue consistently with 95% confidence.

5.2 Selectivity

30 negative control meat samples were evaluated in an independent laboratory and o of these negative samples tested positive with the Premi[®]Test.

5.3 Cross reactivity

Premi[®]Test does cross-react with the following drugs spiked at a concentration of 100 ppb:

| | |
|----------------|--|
| β-lactams: | Amoxicillin, Ampicillin, Cephalirin, Cloxacillin |
| Tetracyclines: | Chlortetracycline, Tetracycline |
| Macrolides: | Erythromycin, Tylosin, Tilmosicin |
| Glycopeptides: | Virginiamycin |
| Sulfonamides: | Sulfadiazine, Sulfadimethoxine, Sulfathiazole |
| Miscellaneous: | Dapsone |

Premi[®]Test does cross-react with the following drugs spiked at a concentration of 50 ppb:

| | |
|----------------|--------------------------------------|
| β-lactams: | Amoxicillin, Ampicillin, Cephalirin, |
| Macrolides: | Erythromycin, Tylosin, Tilmosicin |
| Miscellaneous: | Dapsone |

Premi[®]Test does cross-react with the following drugs spiked at a concentration of 10 ppb:

| | |
|----------------|-------------------------|
| β-lactams: | Amoxicillin, Ampicillin |
| Miscellaneous: | Dapsone |

6. Test principle

- Premi[®]Test is based on the inhibition of the growth of *Bacillus Stearothermophilus var. calidolactis*, a thermophilic microorganism which is very sensitive to many antibiotics and sulphonamide residues. A standardized number of spores is imbedded in an agar medium with selected nutrients and a pH indicator (bromocresol purple).
- When muscle tissue fluid containing no antimicrobials or non-inhibitory substances is added to the Premi[®]Test and heated at 64 °C, the spores will grow, multiply and produce acid. This process is characterized by a visible change in the color of the pH indicator from purple to yellow in the ampoule.
- When muscle tissue fluid containing antimicrobial residues or inhibitory substances are present at or above the detection sensitivity, no growth will occur and there will be no acid production. As a result, there will be no change in the color of the pH indicator and the color of the ampoule will remain purple.

7. Caution

This test is extremely sensitive to antibiotics, disinfectants and other inhibitory substances. Any contamination with such materials should be avoided at all times. It is advised to wash hands thoroughly before starting the test procedure, and to dry the hands with tissue paper or a clean towel. Disposable hygienic hand gloves can also be used.

8. Operating instructions Premi[®]Test

8.1 Before starting the test procedure

- Premi[®]Test: take the number of ampoules needed from storage and let them warm up to room temperature. (18-25 °C / 64 - 77°F)
- Premi[®]Test incubator: Allow the Premi[®]Test incubator or water bath to reach the required temperature of 64°C ± 1.0°C (147°F).
- Samples: Take approximately 4 grams of fresh, lean bovine muscle tissue and use a meat press to extract approximately 250 µl of meat fluid. The juice can be obtained somewhat easier, when heating the meat press first (e.g. by placing it in an incubator of 50°C (122°F) some hours prior testing)
- Number of samples: The number of meat-samples that can be tested in one series is 9. The Premi[®]Test Incubator contains 10- holes. One hole should be reserved for the negative control.
- Wash hands thoroughly before starting the test procedure.

| | |
|--|--|
| | |
| | |

| | |
|--|--|
| | |
| | |

| | |
|--|--|
| | |
| | |

| | |
|--|--|
| | |
| | |

- 8.2 Testing procedure**
- Cut out the required number of ampoules using a pair of scissors. Do not tear off ampoules, as this may damage the foil the adjacent ampoules. Any test with torn foil must be used immediately or discarded.
- Mark the ampoules with a meat-sample identification code.
- Remove foil carefully from the required amount of ampoules (do not open more ampoules than needed: after opening, the test can not be stored any longer).
- Pick up a disposable pipette tip without touching the ends and place it on the syringe.
- Use a new disposable tip on the syringe for each sample.
- Pipette 100 +/- 5 µl of each muscle tissue juice sample to be tested into the labeled ampoules. Depress the plunger of the syringe till the first stop, dip the pipette tip into muscle tissue juice and let the plunger return slowly to its original position. Transfer the meat juice onto the agar of the labeled ampoule by depressing the plunger slowly and completely. Do not disturb the agar.
- Pre-incubate the test 20 minutes at room temperature (18 - 25 °C / 64 - 77 °F)
- Remove the muscle tissue juice by gently emptying the test, by putting the ampoule(s) upside down.
- Wash the test two times, by gently filling and emptying the test ampoule with demineralized (demi) water.
- Remove the last traces of water carefully from the test ampoule, by putting the ampoule up site down on a tissue.
- Close the test ampoule with the foil supplied with the kit to avoid evaporation during incubation.
- Incubate the test ampoule in the Premi[®]Test incubator or water bath (64 °C / 147 °F) until the negative control changes color from purple to yellow. (See Section 6.3)
- When the negative control has changed color (after approximately 3 hours), the ampoules can be read with the Premi[®]Scan.

8.3 Guideline to Determine the Removal of Samples from the Incubator

- Interpreting result for screening by eye (color change) [NOTE: The EVALUATION WAS NOT INTENDED TO COVER THE USE OF THE EYE, but is included as a guideline for triggering the point at which the test samples should be removed from the incubator to be read by the Premi[®]Scan. GUIDELINE FOR TRIGGERING THE POINT AT WHICH THE TEST SAMPLES SHOULD BE REMOVED FROM THE INCUBATOR TO BE READ BY THE PREMISCAN]
- Read the color only from the bottom 2/3 part of the ampoule.
- A partly or completely yellow color of the agar medium indicates the absence of antibiotic residues in the corresponding meat sample at or above the levels indicated above. Such a result should be regarded as a negative result (reported as not found).
- When the color change is "doubtful", incubate 5 minutes longer and read again. When the result is still not clear, test the sample again.

9. Operating Instructions Premi[®]Scan

- Place the Premi[®]Scan in the CD drive; the program will start up automatically.
- Accept the license agreement by clicking on OK and the program will be ready to use.
- Place the Kodak Q-60 Color chart, delivered with the Premi[®]Scan, on the scanner (HP Scanjet 7400 or 8200).
- Calibrate the Premi[®]Scan with the Kodak Q-60 Color chart, by following the instructions of the software.
- Put the ampoules to be tested in the templates, delivered with the Premi[®]Scan.
- Place the black cloth over the templates.
- The Premi[®]Scan recognizes the position of the ampoules automatically and determines the result of the ampoules.
- Results can be printed.

10. Storage and Precautions

- The Premi[®]Test ampoules should be stored in a dry and cool place (preferable below 8 °C / 46 °F). Attention: DO NOT FREEZE! Ampoules should not be used after the expiration date, indicated on the label.
- Anything that may come in contact with the meat (hands, petri-dish etc.) must be cleaned thoroughly, rinsed with tap water and dried with a clean (paper) towel to avoid contamination with inhibiting substances like antibiotics, other drugs, cleaning agents or disinfectants.
- Store the meat below 8 °C (46 °F) if it is not possible to carry out the test immediately. If it is not possible to test the meat within 48 hours, the sample should be stored at or below -18 °C.

11. Technical Assistance

For assistance in use, interpretation, or troubleshooting unanticipated results of this test contact sales@r-biopharm.de. To facilitate the use of Premi[®]Test there is a PowerPoint document available for all users.

12. Limited Liability

Premi[®]Test is a screening test and as such 100% accuracy of the test results cannot be guaranteed. Besides, the assessment of color, in particular that of a yellow/purple result, may differ from person to person. In cases where severe consequences are involved for the user, test results should be confirmed by a validated comprehensive analytical method (like HPLC or LC-MS MS). The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product. The above test has received AOAC approval for this test kit use in North America.

| | |
|--|--|
| | |
| | |

| | | |
|--|---|---|
| R-Biopharm AG Postanschrift / Postal Address: An der neuen Bergstraße 17 64297 Darmstadt, Germany Sitz / Corporate Seat: Pfungstadt Tel.: +49 (0) 61 51 - 81 02-0 Fax: +49 (0) 61 51 - 81 02-40 E-mail: info@r-biopharm.de www.r-biopharm.com | Aufsichtsrat/Supervisory Board: Dr. Ralf M. Dreher (Vorsitzender/Chairman) Vorstand/Board of Management: Christian Dreher (Vorsitzender/Chairman), Dr. Hans Frickel, Jochen Hirsch, Dr. Peter Schubert Handelsregister/ Commercial Register: Amtsgericht Darmstadt, HRB 8321 Sitz der Gesellschaft/ Corporate Seat: Pfungstadt USt.ID-Nr. / VAT-No.: DE 111 657 409 |  |
|--|---|---|