r-biopharm®

RIDA® QUICK Zearalenon RQS

REF R5504

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RIDA®QUICK Zearalenon RQS

Brief information

RIDA®QUICK Zearalenon RQS (Art. No. R5504) is an immuno-chromatographic test for the quantitative determination of zearalenone in corn (see chapter 1. Intended use).

The test kit contains 20 test strips, each of which can be used for one analysis. All reagents required for the assay are contained in the test kit. Results are evaluated with the RIDA®SMART APP software (Art. No. ZRSAM1000) or the RIDA®QUICK SCAN reader.

Sample preparation: homogenization and extraction

Time requirement: sample preparation (for 10 samples)... approx. 10 min
test implementation (incubation time)..... 5 min

Limit of detection: approx. 75 µg/kg (ppb)
(evaluation with RIDA®QUICK SCAN)
approx. 50 µg/kg
(evaluation with RIDA®SMART APP)

Measuring range: 50 - 1000 µg/kg / 75 - 500 µg/kg
Quantitative evaluation with the RIDA®SMART APP/
RIDA®QUICK SCAN

Specificity: The RIDA®QUICK Zearalenon RQS test reacts with zearalenone in corn samples.

Further information is contained in the validation report.

The specificity of the RIDA®QUICK Zearalenon RQS test was determined by analyzing the cross reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

Related products and accessories for zearalenone determination

RIDASCREEN® Zearalenon (Art. No. R1401)

RIDASCREEN®FAST Zearalenon (Art. No. R5502)

RIDASCREEN®FAST Zearalenon SC (Art. No. R5505)

RIDA®SMART APP software (Art. No. ZRSAM1000)

1. Intended use

RIDA®QUICK Zearalenon RQS (Art. No. R5504) is a quantitative immuno-chromatographic test in strip format for the quantitative determination of zearalenone in corn.

2. General information

The mycotoxin zearalenone is formed by fungi of the genus *Fusarium*. Zearalenone is a phytohormone which displays, apart from its anabolic properties, mainly estrogenic effects. Because of its estrogenic properties, zearalenone may induce fertility disorders in animals with clinical signs of hyperestrogenism - an aspect of a disease which although reported mainly in hogs, is described in other species such as cow, horse and sheep, too.

The potential health risk for humans induced by this mycotoxin, which is taken up with foods of vegetable or animal origin, is extensively discussed.

3. Test principle

The immunochromatographic test in the form of a test strip is based on an antigen-antibody reaction. A specific anti-zearalenone antibody detects zearalenone in the sample. During incubation of the test strip, a band pattern (test line / control line) forms that is used to determine the concentration of zearalenone. The intensity of the test line depends on the zearalenone concentration of the sample. It increases as the zearalenone concentration increases. It must be possible to detect the control line after the reaction has taken place in order to verify the functioning of the test. The control line weakens as the zearalenone concentration in the sample increases.

The test strip is evaluated using the RIDA®SMART APP software (Art. No. ZRSAM1000) or the RIDA®QUICK SCAN reader.


Fig. 1: RIDA®QUICK Zearalenon RQS test strip.

4. Reagents provided

Each kit contains sufficient materials for max. 20 determinations. Each test kit contains:

Component	Cap color	Format	Volume
20 x Test strip	-	Ready to use	One determination per test strip (separately packed)
1 x Mobile solvent	White	Ready to use	13 mL
1 x RIDA®SMART APP cover	-	Ready to use	Lot-specific evaluation cover for use with RIDA®SMART APP

5. Reagents required but not provided

5.1 Equipment

- Gloves
- Laboratory grinder
- Balance
- Graduated cylinder
- Optional: Shaker (e.g. Rock-it 360, Trilogy® Analytical Laboratory Inc. or equivalent)
- Optional: Centrifuge ($\geq 2,000$ g)
- 100 μ L and 600 μ L micropipettes
- Latest version of the RIDA®SMART APP software (Art. No. ZRSAM1000) activated and installed on a compatible smartphone or RIDA®QUICK SCAN

5.2 Reagents

- Distilled water (dist. Water) or deionized water
- 70 % methanol (e.g. 70 mL 100 % methanol + 30 mL distilled water), exact adjustment necessary

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

7. Storage instructions

Store the kit at 2 - 8 °C (36 - 46 °F). Do not freeze the test strips or any test kit components.

The test strips are sensitive to moisture. Moist test strips can have a negative impact on the test result. For this reason, it is important to protect the strips against moisture. Do not remove the test strips from the test strip packaging until immediately prior to use in the test.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Sample preparation

Bring all components necessary for sample preparation to room temperature (20 - 25 °C / 68 - 77 °F) before use in the test and perform sample preparation at room temperature. Do not allow methanol to evaporate.

The sample should be stored in a cool place, protected against light.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

- Weigh 10 g of ground and mixed sample into a suitable container and add 20 mL 70 % methanol.
- Close the tube and shake the sample vigorously approx. 5 min (manually or with shaker / vortex).

- Let the solution settle for 3 min or centrifuge.
- Mix 600 µL of the temperate mobile solvent with 100 µL of clear supernatant.
- Thereof use 100 µL in the test.

Note

Please avoid methanol evaporation during sample preparation.

If necessary, the sample weight can be increased, e.g. 20 g sample in 40 mL 70 % methanol. In this case the volume of methanol has to be adapted.

9. Test procedure

Bring all components necessary for sample preparation to room temperature (20 - 25 °C / 68 - 77 °F) before use in the test and perform sample preparation at room temperature.

- Remove a test strip from its foil pouch.
- Apply 100 µL of the sample solution (see chapter 8. Sample preparation) on the application area of the test strip (see Fig. 1).
- Incubate the test strip at room temperature.
- Evaluate the result after 5 min (+ max. 15 seconds) incubation with the RIDA®SMART APP or RIDA®QUICK SCAN reader.

10. Evaluation

The left line in the reaction area is the test band (test line, see fig. 1). Its appearance and intensity depends on the zearalenone concentration of the sample. The sample is contaminated with zearalenone, if the control band (control line) is visible and the test band (test line) is also visible.

The right line in the reaction area is a control band (control line) and must be present after each test procedure. If this line is not visible, the test was not performed correctly or the reagents were not in order. In this case, the test should be repeated with a new test strip. If the right band is missing again, please inform your local distributor.

Please note: If the sample is highly contaminated with zearalenone, the control line will be only faintly visible. In this case the test result might be reported as “invalid result“ (RIDA®SMART APP) or “invalid“ (RIDA®QUICK SCAN reader).

To obtain comparable results, the test strips must always be evaluated after 5 min.

10.1 Evaluation with RIDA®SMART APP

First, please read the instructions for use of the RIDA®SMART APP software (Art. No. ZRSAM1000) attentively. The test strip evaluation has to be performed by using the RIDA®SMART APP software according to these guidelines. Please find a short description (Quick Starter Guide) of how to use the RIDA®SMART APP on our website www.r-biopharm.de.

For the test strip evaluation the RIDA®SMART APP cover is needed which is included in the test kit. The RIDA®SMART APP cover contains lot-specific information and must be used only with the appropriate lot of the test kit. For further information please check the Quality Assurance Certificate enclosed in the kit. Please make sure the correct RIDA®SMART APP evaluation cover is present prior to starting the test procedure and test strip evaluation with the RIDA®SMART APP software.

RIDA®SMART APP software applications for evaluation:

Matrix	Measuring range	RIDA®SMART APP application
Corn	50 - 1000 µg/kg	50-1000ppb

10.2 Evaluation with RIDA®QUICK SCAN

First, carefully read the user guide for the RIDA®QUICK SCAN and proceed as instructed. Point 3. contains the description for evaluating the test strips.

To obtain optimal results every test kit contains lot-specific standard curve parameters. These parameters can be found on the Quality Assurance Certificate and have to be used only with this specific lot. The Quality Assurance Certificate contains parameters for several standard curves. Depending on the extraction method used, select the corresponding parameters.

To enter the parameters manually into the RIDA®QUICK SCAN reader use the plain text information on the Quality Assurance Certificate. To enter the data with the barcode scanner use the 2D barcodes. In both cases, the checksum shown on the RIDA®QUICK SCAN screen needs to be identical to the respective checksum on the certificate (number below the 2D barcode).

RIDA®QUICK SCAN evaluation methods:

Measuring range	RIDA®QUICK SCAN method
75 - 500 µg/kg	Zea QUANT or Zea75-500ppb

11. Sensitivity

This test determines zearalenone contaminations in the measurement range of 50 - 1000 µg/kg in combination with the RIDA®SMART APP software and in the measurement range of 75 - 500 µg/kg in combination with the RIDA®QUICK SCAN reader.

12. Further application notes

Further application notes are available on request. Please contact your local distributor.

For further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2011-10-24	Release version
2016-06-30	General revision
2021-12-07	Current version General revision Changes made: – 3. Test principle: Figure from previous chapter 10. Evaluation listed and updated here (new test strip housing)

Explanation of symbols

General symbols:

	Follow the instructions for use
	Batch number
	Expiry date (YYYY-MM)
	Storage temperature
	Article number
	Number of test determinations
	Manufacturing date (YYYY-MM)
	Manufacturer + address

Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.