RIDASCREEN® Peanut

REF R6811

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RIDASCREEN® Peanut

Brief information

RIDASCREEN® Peanut (Art. No. 6811) is a sandwich enzyme immunoassay for the quantitative analysis of peanut and peanut protein in food validated for the method (see chapter 1. Intended use).

All reagents required for the enzyme immunoassay, including standards, are contained in the test kit. The test kit is sufficient for a maximum of 96 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

The RIDASCREEN® Peanut is certified at AOAC-RI (PTM No. 112102).

- Sample preparation: homogenization, extraction and centrifugation
- Time requirement: sample preparation (for 10 samples)... approx. 20 min
test implementation (incubation time)..... 50 min
- Standard material: NIST Reference Material 2387 (peanut butter)
(Protein content: 22.2 ± 0.5 g / 100 g)
- Limit of detection: 0.15 mg/kg (ppm) peanut
(corresponds to 0.033 mg/kg (ppm) peanut protein)
for further information on the statistical evaluation see
the validation report
- Limit of quantification: 0.75 mg/kg (ppm) peanut
(corresponds to 0.166 mg/kg (ppm) peanut protein)
- Specificity: The antibodies used specifically detect peanut proteins. The system shows no cross reactivity to 91 tested food commodities.

Further information are contained in the validation report.

Cross reactivities of the antibodies used for this test kit have been determined for the pure food (e.g. corn flour). In composed / processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spiking experiments (see chapter 13. Limits of the method).

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our Good ELISA Practice brochure. It lists minimum standards and conditions that are required when using test kits of R-Biopharm AG to perform ELISA analysis. The brochure can be retrieved, printed and downloaded from the website

<https://food.r-biopharm.com/media/technical-guides/>.

Related products for peanut determination

Bioavid lateral flow Erdnuss / Peanut (Art. No. BL606-10; BL606-25)

Bioavid lateral flow Erdnuss / Peanut incl. hook line (Art. No. BLH706-20)

SureFood® ALLERGEN Peanut (Art. No. S3603)

SureFood® ALLERGEN 4plex Peanut / Hazelnut / Walnut + IAC (Art. No. S3402)

1. Intended use

RIDASCREEN® Peanut (Art. No. R6811) is a sandwich enzyme immunoassay test developed for the quantitative analysis of peanut and peanut protein in food. Due to the large number of different food products, the following samples were examined as representative for different food product categories within the scope of the test development: puffed rice cereals, granola bar, cookies, ice cream, trail mixes and milk chocolate. It can be assumed that the assay is also suitable for the analysis of other foods, however, this must be verified by the user before applying the test kit to that food category.

For detailed results and further information on validation data with other food matrices, please refer to the validation report.

2. General information

Although peanut (*Arachis hypogaea*) botanically belongs to the legume family and not to the nuts family, its taste and usual consumption is similar to that of tree nuts. Peanut is, so to speak, one of the most popular "nuts" consumed worldwide. They are very common ingredients in chocolate products and are also consumed as snacks and in baked goods such as cookies.

The peanut protein content of 25 % in peanuts is very high. Many of these proteins are known to be allergenic, such as the arachins and the conarachins, which are present in relatively high amounts. Peanut is also one of the eight foods that are responsible for the most common cause of food allergies. Allergies to peanuts are among the most severe allergies and can lead to life-

threatening anaphylactic shocks. The consumption of less than one milligram of peanut can cause allergic reactions in very sensitized people. Therefore, labeling of peanut-containing products is mandatory in the USA and the European Union as well as in many other countries. The main danger for sensitized persons is the contamination of peanuts in food products that are supposed to be peanut-free (e.g. peanut traces in chocolate that should not contain peanuts). It is therefore important to test food that may have been contaminated during the production process for the presence of unintentional peanut traces.

In 2019, Taylor et al. (2019 VITAL Scientific Expert Panel Recommendations, The Allergen Bureau Limited (<http://allergenbureau.net>)) re-evaluated the VITAL values for peanut and found that the ED01 value - the dose at which 99 % of peanut-allergic individuals will not react - is 0.2 mg of peanut protein. The ED05 value - the dose at which 95 % of peanut-allergic individuals will not react - is 2.1 mg of peanut protein.

3. Test principle

The principle of the test is the antigen-antibody reaction. The wells of the microtiter strips are coated with specific antibodies against peanut proteins. By adding the standard or sample solution to the wells, peanut proteins present in the sample will bind to the specific capture antibodies resulting in the formation of an antibody-antigen-complex. Components not bound by the antibodies are then removed in a washing step. Following the washing step, a solution containing antibody conjugated to peroxidase is added. This conjugate is bound to the Ab-Ag-complex and an antibody-antigen-antibody (sandwich) complex is formed. Any unbound conjugate is then removed in another washing step. A substrate/chromogen solution is added to the wells and incubated. Bound conjugate converts the chromogen into a blue product. A stop solution is added which results in a color change from blue to yellow. The absorbance of the solution which is proportional to the peanut protein concentration in the sample is measured photometrically at 450 nm and expressed as mg/kg peanut or peanut protein.

4. Reagents provided

Each kit contains sufficient materials for a maximum of 96 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		96 wells
Allergen extraction buffer	Green	Concentrate	10x	100 ml
Standard 1*	Transparent	Ready to use	0 mg/kg	1.3 ml
Standard 2*	Transparent	Ready to use	0.75 mg/kg	1.3 ml
Standard 3*	Transparent	Ready to use	1.5 mg/kg	1.3 ml
Standard 4*	Transparent	Ready to use	3.0 mg/kg	1.3 ml
Standard 5*	Transparent	Ready to use	6.0 mg/kg	1.3 ml
Wash buffer	Brown	Concentrate	10x	100 ml
Conjugate	Red	Ready to use		11 ml
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use		13 ml
Stop solution	Yellow	Ready to use		14 ml

* The dilution factor 20, which results after sample preparation, has already been considered for the standard concentrations. Therefore, the peanut concentration of the samples can directly be read from the standard curve.

5. Reagents required but not provided

5.1. Equipment

- Gloves
- Scale (measurement range at least 50 g and precision of ± 0.01 g)
- Laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Centrifuge (at least 2,500 x g) + centrifugal vials with cap (e.g. 50 ml centrifuge tubes from Greiner, Art. No. 227261)
- Shaker
- Water bath (60 °C / 14 °F; for fluctuation range please refer to the instructions of the water bath manufacturer)
- Fluted filter (pore size 8 - 12 μ m)
- Graduated pipettes
- Graduated cylinders
- Variable 20 - 200 μ l and 200 - 1000 μ l micropipettes
- If necessary: a further microtiter plate (e.g. universal binding, breakable MTP from Thermo Fisher Scientific, Art. No. 95029390 or low binding Greiner bio-one, Art. No. 655901)
- If necessary: 8-channel pipette for 100 μ l
- Microtiter plate spectrophotometer (450 nm)

- Optional: RIDASOFT® Win.NET (Art. No. Z9996FF)

5.2. Reagents

- Distilled water (dist. water) or deionized water
- Skim milk powder (SMP; food quality; peanut-free)

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory personnel. The instructions for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

Do not reuse wells of the microtiter strips (coated microtiter plate and pre-plate, if necessary, see chapter 10.2. Test procedure)). Use separate pipette tips for each standard and each sample extract to avoid cross contamination.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

To avoid moisture inside the wells, open the foil bag for withdrawal of microwells only after having reached room temperature (20 - 25 °C / 68 - 77 °F).

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive. Therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- Bluish coloration of the reddish substrate/chromogen prior to test implementation
- Value of less than 1.2 absorbance units ($A_{450\text{ nm}} < 1.2$) for standard 5

9. Sample preparation

Wear gloves before starting and during the assay. Airborne allergens and dirty laboratory equipment lead to contamination of the assay. Therefore, please notice the following recommendations:

- Clean surfaces, glass vials, mincers and other equipment before and after each sample preparation.
- Carry out the sample preparation in a room isolated from the ELISA procedure.

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The Allergen extraction buffer is provided as a 10-fold concentrate and must be diluted prior use. Before dilution of the buffer concentrate, dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. After that, dilute the heated buffer concentrate 1:10 (1+9) with dist. water (e.g. 100 ml buffer concentrate + 900 ml dist. water). The **diluted Allergen extraction buffer (AEB)** is stable at 20 - 25 °C (68 - 77 °F) for approx. 4 weeks or at 2 - 8 °C (35 - 46 °F) for 12 weeks.

Make sure to heat the AEB in time in the 60 °C (140 °F) water bath.

In the following chapter, the following abbreviation is used:

- AEB: final diluted Allergen extraction buffer

9.1. Extraction with AEB and skim milk powder (SMP)

Heat the AEB to 60 °C (140 °F) before sample extraction.

- Homogenize (grind thoroughly to powder and mix well or mix well a solution respectively) the laboratory sample to ensure that a representative test portion can be obtained.
- Transfer 1 g (or in case of liquid samples 1 ml) of homogenized sample to a new vial and add 1 g SMP.
- Add 20 ml (or in case of liquid samples 19 ml) pre-heated AEB (see chapter 9.).

- Mix vigorously (e.g. Vortexer).
- Extract for 10 min at 60 °C (140 °F) in a water bath by shaking occasionally.
- Let the sample cool down in ice water for 3 - 5 min.
- Filter or centrifuge sample for 10 min at > 2,500 x g; if possible at 4 °C (39 °F).
(Alternatively: transfer 2 ml of the extract into a reaction vial and centrifuge at high speed (> 10,000 x g) for 10 min in a microcentrifuge).
- Transfer the supernatant to a new vial.
- If the supernatant is not free of particles after centrifugation, filter the extract.
- The extract (supernatant of centrifugation step or filtrate) should be tested by ELISA on the same day if stored at room temperature (20 - 25 °C / 68 - 77 °F). It can be stored for up to two days in a well-closed container at 2 - 8 °C. Unused extracts can be stored at - 20 °C (- 4 °F) for one month.

Remark:

As an alternative to weighing the SMP to the sample, the SMP can be added to the AEB (1 g SMP per 20 ml AEB). The shelf life of the SMP-AEB is only 1 day.

In case of samples with concentration above the calibration range, further dilutions should be made with the SMP-AEB (1 g SMP per 20 ml AEB).

10. Test implementation

10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **wash buffer** is provided as a 10-fold concentrate. Before use, the buffer has to be diluted 1:10 (1+9) with dist. water (i.e. 100 ml buffer concentrate + 900 ml dist. water). Prior to dilution, dissolve eventually formed crystals by incubating the buffer in a water bath at 37 °C (99 °F). The diluted buffer is stable at 20 - 25 °C (68 - 77 °F) for 4 weeks.

Components should be stored immediately at 2 - 8 °C (35 - 46 °F) when no longer required.

10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than three strips (24 wells) at a time. If more than three strips are needed, a second uncoated plate (see chapter 5.1.) should be used as a

pre-plate to avoid a time shift over the microtiter plate. All standards and samples are pipetted into the uncoated plate (at least 150 µl per well) and then exactly 100 µl are quickly transferred to the coated microtiter plate with an 8-channel pipette.

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 100 µl of each standard or sample (prepared according to chapter 9.) in duplicate to the wells and incubate for 20 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl diluted wash buffer (see chapter 10.1.) and pour out the liquid as before. Repeat three more times (a total of four wash cycles).
4. Add 100 µl of the conjugate to each well and incubate for 20 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times in a row) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl diluted wash buffer (see chapter 10.1.) and pour out the liquid as before. Repeat three more times (a total of four wash cycles).
6. Add 100 µl of the reddish substrate/chromogen to each well and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 30 min¹ after addition of stop solution.

¹ Current data show that measurement is possible up to 30 min after addition of the stop solution (not part of the AOAC approval). For the AOAC approval, a measurement within 10 min was specified.

11. Evaluation

A special software, the **RIDASOFT® Win.NET (Art. No. Z9996FF)**, is optional available for evaluation of the RIDASCREEN® enzyme immunoassays. The evaluation should be done using the 4-parameter function.

For the evaluation it should be confirmed, that the quality criteria are fulfilled for the current test run. The shape of the standard curve should closely match the shape of the standard curve printed in the analysis certificate.

The assay can be also evaluated when running in single well per extract. This has no influence on the function of the test kit. A special assay evaluation must be written in the RIDASOFT® Win.NET software for this purpose. It is not present by default. Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. However, it is not consistent with standards like EN 15633-1 and EN 15842 and also does not correspond to the procedure for the AOAC Approval. It should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, a higher result variation will occur when pipetting in single wells per extract.

12. Result interpretation

The test result is calculated in mg peanut per kg food. A conversion factor of 0.222 determined from protein content of 22.2 % of the peanut reference standard material can be used to calculate the peanut protein concentration in mg/kg. Multiply the result obtained in mg peanut per kg food by 0.222 to obtain the peanut protein content in the food matrix.

Results between LOD and LOQ indicate a low peanut concentration in the sample. Calculated results show a high uncertainty in this area due to the method's high variation below LOQ. Therefore, such results should not be reported with a quantitative value, but qualitative as "< LOQ".

A result below the LOD does not exclude an allergen contamination below the detection limit of the assay, or that other peanut components, such as lipids, may be present in a sample. The result should be reported accordingly.

A further dilution and new detection of samples is recommended for absorbance values ($A_{450\text{ nm}}$) > standard 5. In case of a further dilution, the additional dilution factor must be taken into account when calculating the peanut concentration. Further dilutions should be made with the SMP-AEB (1 g SMP per 20 ml AEB).

Compared to the certificate, higher absorbance values ($A_{450\text{ nm}}$) for the standard curve, especially for the zero standard, may be a result of insufficient washing or peanut contamination.

13. Limits of the method

Test results may vary depending on the sample matrix, the actual test procedure and the laboratory environment.

Detection and quantification limits depend on the respective sample matrix, the degree of processing and the extraction method.

The protein content and protein composition may vary in different peanut varieties. Therefore, different peanut varieties may give different results, as the test was calibrated against exemplary peanut varieties used with the standard material.

Technical limits of the test method are approached outside the designated measurement range resulting in higher variation. This may cause a switch of results between the different areas of the calibration curve especially at the test characteristic boundaries (LOD, LOQ, upper limit of measurement range).

An incorrect weight of the sample to be analyzed will have a 1:1 effect on the measurement result (e.g. a 10 % higher concentration is measured with a weigh in of +10 %). A sufficient accuracy is given with a fluctuation of max ± 1 %.

For detailed results and further information for other food matrices, please refer to the current validation report. In addition, data on individual foods may be available from comparative laboratory tests and inter-laboratory comparisons.

For the present ELISA, only individual, exemplary foods from different product categories could be validated due to the large number of foods. When analyzing a non-validated matrix, it is recommended to verify the results obtained by means of spike experiments. If necessary, a validation of the sample matrix of interest will need to be performed.

Due to the multitude of food types, matrix effects cannot be excluded. These can lead to false-positive / increased results, but also reduce or suppress a correct reaction. Such matrix effects are independent of the specificity of the antibody used in the test and can be made visible by spiking experiments.

In processed foods (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented and this may have an impact on the recovery and assay results.

Cross reactivities are side reactions of the antibody used for preparing the test kit with antigen showing similar epitopes as the investigated analyte. These appear especially with antigens from closely related species. In contrast to matrix effects, it is a specific reaction of the antigen with the antibody. The

antigen structures are subject to similar influences (e.g. by heating or drying) as the actual analyte. Therefore, cross reactivities may also appear after food processing in single case or are lost.

For evaluation of the cross reactivity only one representative sample was analyzed, other samples may show a different result. All analyzed cross reactivities are described in the validation report.

14. Recommendation

In order to ensure a high analytical performance we recommend:

- To comply with the general quality assurance requirements for laboratories as listed in standards like EN 15633-1 and EN 15842 (e.g. performing duplicate determinations).
- To pre-flush pipette tips with standard or sample extract prior to pipetting.
- To carry along test controls for quality control and to ensure an accurate and correct test procedure. Peanut-free and peanut-containing (naturally contaminated or spiked) samples should be used. An example of a spike experiment is given in the validation report.
- In case of extremely acidic or basic samples, adjust the sample's pH value (pH 6.5 - 7.5) to neutral prior to extraction.
- To perform a PCR (e.g. SureFood[®]) for confirmation of the result.
- To contact sales@r-biopharm.de if automates (e.g. ThunderBolt[®] / Bolt[™]) are used.

Further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2021-02-01	Release version
2021-12-16	Current version Changes made: <ul style="list-style-type: none">– New statistical evaluation of LoD (0.15 mg/kg) according to AOAC– General linguistic revision– Adaptations in chapter 1. Intended use– Modification of the disposal clause in chapter 6.– Chapter 10.2. Test procedure: time span adding stop solution until measuring increased from 10 min to 30 min– Addition in chapter 11. Evaluation (moved from chapter Recommendation)– Revision of chapter 11. Evaluation, 12. Result interpretation, 13. Limits of the method, 14. Recommendation and 15. Further application notes

Explanation of symbols

General symbols:

	Follow the instructions for use
	Batch number
	Expiry date (YYYY-MM)
	Storage temperature
	Article number
	Number of test determinations
	Manufacturing date (YYYY-MM)
	Manufacturer + address

Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.