**CONGEN**

# **SureFood<sup>®</sup> ALLERGEN Gluten**

Art. No. S3606  
100 rxn

## **User Manual**

**July 2022**

 **Inhalt**

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## **Content**

1	General Information .....	11
1.1	Description .....	11
1.2	Limit of Detection and Limit of Quantification .....	11
1.3	DNA-preparation .....	12
1.4	Kit components and storage .....	12
1.5	Additionally required equipment and materials .....	12
1.6	Setup .....	12
1.7	Detection channel Set-up .....	13
2	Qualitative Analysis .....	14
2.1	Protocol .....	14
2.1.1	Preparation of the master-mix .....	14
2.1.2	Preparation of the real-time PCR-mix .....	14
2.2	Interpretation of results .....	15
3	Quantitative Analysis .....	16
3.1	Protocol .....	16
3.1.1	Preparation of the master-mix .....	16
3.1.2	Preparation of the standard DNA dilutions .....	16
3.1.3	Preparation of the real-time PCR-mix .....	17
3.2	Interpretation of results .....	17
4	Further Information .....	18
4.1	Product Information .....	18
4.2	Technical Support .....	18
4.3	Distribution and Ordering .....	18

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## **1 General Information**

### **1.1 Description**

The SureFood® ALLERGEN Gluten is a real-time PCR for the direct, qualitative and / or quantitative detection of specific gluten-containing cereals (wheat such as spelt and khorasan wheat, rye, barley, oats) DNA sequences according to directive (EC) 1169/2011.

Each reaction contains an internal amplification control (IAC). If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). In addition spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

For the quantitative determination the use of the laboratory reference material SureFood® QUANTARD Allergen 40 (Art. No. S3301) containing 40 mg wheat / kg food sample and a detection system including a standard series is required.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of two fluorescence emissions at the channels FAM and VIC/HEX at the same time. The technical verification of instruments was performed on Roche LightCycler® 2.0, Roche LightCycler® 480 II, Roche cobas® z 480 Analyzer, Qiagen Rotor-Gene Q, Applied Biosystems 7500, Bio-Rad CFX96, Bio Molecular Systems MIC, LTF MyGo Pro, Agilent AriaDx and Agilent Mx3005P.

In general this procedure is applicable for all homogenous food samples. The quantitative determination is inapplicable for swabs, sponges and inhomogeneous fluids/liquids.

The method is validated for quantitative determination of 1 mg to 400 mg allergenic substance / kg food sample using the reference material SureFood® QUANTARD Allergen 40 containing 40 mg allergenic substance / kg food sample.

### **1.2 Limit of Detection and Limit of Quantification**

The SureFood® ALLERGEN Gluten PCR has a limit of detection of  $\leq 0.4$  mg / kg using SureFood® PREP Advanced, protocol 1.

The limit of quantification is 1 mg / kg using SureFood® PREP Advanced, protocol 1.

The limit of detection and quantification were determined using the SureFood® QUANTARD Allergen 40 (matrix: corn flour). The limit of detection and quantification of the complete method (DNA extraction and real-time PCR) depends on sample matrix, processing grade, DNA-preparation and DNA-content.

The SureFood® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

**1.4 DNA-preparation**

For DNA-preparation the use of SureFood® PREP Advanced (Art. No. S1053), protocol 1 or SureFast® Mag PREP Food (Art. No. F1060) is recommended.

**1.5 Kit components and storage**

<b>Kit Code</b>	<b>Reagent</b>	<b>Amount</b>	<b>Lid Color</b>
1	Reaction Mix	2 x 1100 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue
4	Dilution Buffer	1 x 1800 µl	White
5	Standard DNA	1 x 200 µl	Dark Blue

**Store all reagents at –20°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.**

**Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.**

**1.6 Additionally required equipment and materials**

- DNA-Extraction kit (e.g. SureFood® PREP Advanced Art. No. S1053 or SureFast® Mag PREP Food Art. No. F1060)
- reference material for quantification (SureFood® QUANTARD Allergen 40, Art. No.: S3301)
- real-time PCR instrument with two detection channels (510 nm and 580 nm)
- real-time PCR consumables (plates, tubes, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

**1.7 Setup**

	<b>Blockcycler &amp; R-Biopharm RIDA®CYCLER</b>	<b>Rotorcycler &amp; LTF MyGo Pro</b>
Initial Denaturation (HOLD) Cycles	5 min, 95°C 45	1 min, 95°C 45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

**1.8 Detection channel Set-up**

<b>Real-time PCR device</b>	<b>Detection</b>	<b>Detection channel</b>	<b>Quencher</b>	<b>Note</b>
<b>Agilent Mx3005P</b>	gluten-containing cereals	FAM	+	
	IAC	HEX	+	
<b>Agilent AriaMx /Dx</b>	gluten-containing cereals	FAM	+	
	IAC	HEX	+	
<b>Applied Biosystems 7500</b>	gluten-containing cereals	FAM	None	Check the passive reference option ROX is none.
	IAC	VIC	None	
<b>Bio-Rad CFX96/ Dx</b>	gluten-containing cereals	FAM	+	
	IAC	VIC/HEX	+	
<b>LTF MyGo Pro</b>	gluten-containing cereals	FAM	+	
	IAC	VIC	+	
<b>Qiagen Rotor-Gene Q</b>	gluten-containing cereals	green	+	<b>Note:</b> Please use only 0.1 ml reaction tubes. The gain settings must be set to 5 (factory default) for all channels.
	IAC	yellow	+	
<b>R-Biopharm RIDA®CYCLER</b>	gluten-containing cereals	green	+	
	IAC	yellow	+	
<b>Roche LightCycler®2.0</b>	gluten-containing cereals	530	None	The SureCC Color Compensation Kit II (Art. No. F4010) is required. Check the Seek Temperature is 58°C.
	IAC	560	None	
<b>Roche LightCycler® 480 II</b>	gluten-containing cereals	465-510	+	
	IAC	533-580	+	
<b>Roche cobas® z 480 Analyzer</b>	gluten-containing cereals	465-510	+	
	IAC	540-580	+	

## 2 Qualitative Analysis

### 2.1 Protocol

#### 2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay.

Recommended control reactions for the specific PCR assay: negative control, extraction control, positive control. The reaction mix contains an internal amplification control (IAC) per reaction.

#### Reactions needed for the qualitative detection of gluten-containing cereals:

3 reactions for controls (1x no-template control, 1x extraction control, 1 positive control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

#### Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

**Mix each master-mix well and centrifuge shortly before use.**

#### 2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

## **2.2 Interpretation of results**

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

A sample is stated **positive**, if the sample DNA shows amplification in the detection system or the obtained Cp value is  $\leq 35$ . High amplicon concentrations can result in a weak or absent signal of the internal amplification control (IAC).

A Cp value for the internal amplification control (IAC) is not needed to obtain a positive result of the Positive control / Standard DNA.

A sample is stated **negative**, if the sample DNA shows no amplification in the detection system or the obtained Cp value is  $> 35$ . The internal amplification control (inhibition control) of the sample has to be **positive** (VIC/HEX) with a shift in Cp-Value  $\leq 2$  compared to the negative control. If the sample DNA in the VIC/HEX-Channel shows **no amplification** or a shift in Cp-value  $> 2$  compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific gluten-containing cereals PCR assay.

Alternatively, a self-selected reference sample at the detection limit can be used to assess the sample. This reference sample is used outside CONGEN's quality and evaluation guidelines.

### **In general:**

Samples tested negative still could contain an allergen contamination below the limit of detection of the assay, or they might contain other allergen components like proteins or lipids for example.

### 3 Quantitative Analysis

#### 3.1 Protocol

##### 3.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay as well as for the inhibition control. Recommended control reactions for the specific PCR assay: negative control, extraction control, standard curve as positive control. The reaction mix contains an internal amplification control (IAC) per reaction.

##### Reactions needed for the quantitative detection of gluten-containing cereals:

5 reactions for the standard curve

2 reactions for controls (1x no-template control, 1x extraction control)

2 reactions for reference material (SureFood® QUANTARD Allergen 40 DNA)

For each sample: at least 1 reaction for each sample DNA

##### Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

Mix each master-mix well and centrifuge shortly before use.

##### 3.1.2 Preparation of the standard DNA dilutions

Use the undiluted Standard DNA (**Code 5**) and dilute the Standard DNA (**Code 5**) in 1:10 steps in Dilution Buffer (**Code 4**) in order to prepare different concentrations. Prepare a dilution series of 4 steps. Prepare 4 reaction tubes (labelled S2 to S5) and add 45 µl Dilution Buffer each.

The following procedure is recommended:

Standard	Dilutions	Copy number per µl	Final copy number per reaction*
S1	Standard DNA undiluted	100,000 copies	500,000 copies
S2	45 µl Dilution Buffer + 5 µl DNA of S1	10,000 copies	50,000 copies
S3	45 µl Dilution Buffer + 5 µl DNA of S2	1,000 copies	5,000 copies
S4	45 µl Dilution Buffer + 5 µl DNA of S3	100 copies	500 copies
S5	45 µl Dilution Buffer + 5 µl DNA of S3	10 copies	50 copies

\*Note: 5 µl of standard DNA are used for each calibration point. The final copy number per reaction is to be entered in the analysis software of the real-time PCR detection system.

### **3.1.3 Preparation of the real-time PCR-mix**

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of the SureFood® QUANTARD 40 Allergen DNA and the standard dilutions into the designated tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

### **3.2 Interpretation of results**

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

A sample is stated **positive and can be quantified**, if the sample DNA shows amplification in the detection system or the obtained Cp value is  $\leq 35$ . High amplicon concentrations can result in a weak or absent signal of the internal amplification control (IAC).

A Cp value for the internal amplification control (IAC) is not needed to obtain a positive result of the Positive control / Standard DNA.

A sample is stated **negative**, if the sample DNA shows no amplification in the detection system or the obtained Cp value is  $> 35$ . The internal amplification control (inhibition control) of the sample has to be **positive** (VIC/HEX) with a shift in Cp-Value  $\leq 2$  compared to the negative control. If the sample DNA in the VIC/HEX-Channel shows **no amplification** or a shift in Cp-value  $> 2$  compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific gluten-containing cereals PCR assay.

The calculation of mg allergenic substance/kg food sample can be done for samples showing an amplification curve in the detection system with Cp values  $< 35$ . Mark the standards, controls and samples and make the evaluation according to the analysis program recommended by the real-time PCR instrument manufacturer. The value for the slope of the standard curve has to be between -3.1 and -3.6 and the correlation coefficient  $R^2 > 0.98$ . In case of different values for the standard curve, it should not be used for calculation.

#### **In general:**

Samples tested negative still could contain an allergen contamination below the limit of detection of the assay, or they might contain other allergen components like proteins or lipids for example.

**For a final calculation of mg allergenic substance / kg food sample following formula is introduced:**

$$X[mg/kg] = \frac{10^{\left(\frac{Cp_{sample}}{s}\right)}}{10^{\left(\frac{Cp_{Quantard}}{s}\right)}} * C_{Quantard}$$

s = slope | c = concentration | Cp<sub>Quantard</sub> = Ø of double assay

**Example:**

$$X[mg/kg] = \frac{10^{\left(\frac{28.5}{-3.251}\right)}}{10^{\left(\frac{29.2}{-3.251}\right)}} * 40 = 65.7 mg/kg$$

Cp<sub>sample</sub> = 28.5 | Cp<sub>Quantard</sub> = Ø 29.2 | s = - 3.251 | c = 40 mg/kg

For this example an amount of **65.7 mg gluten-containing cereals / kg food sample** is calculated.

## 4 Further Information

### 4.1 Product Information

- Microsoft Excel template of calculation and detailed information about setup of several real-time PCR devices (Download: [www.congen.de/en/company/downloads](http://www.congen.de/en/company/downloads))
- Validation Report upon request

### 4.2 Technical Support

For further questions please contact your distributor or send an e-mail to [sales@r-biopharm.de](mailto:sales@r-biopharm.de).

### 4.3 Distribution and Ordering

R-Biopharm AG  
An der neuen Bergstrasse 17,  
64297 Darmstadt, Germany  
Phone: +49 (0) 61 51 - 81 02-0  
Fax: +49 (0) 61 51 - 81 02-20  
E-Mail: [orders@r-biopharm.de](mailto:orders@r-biopharm.de)  
[www.r-biopharm.com](http://www.r-biopharm.com)