RIDASCREEN® SET A,B,C,D,E

REF R4101

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RIDASCREEN® SET A,B,C,D,E

Brief information

RIDASCREEN® SET A,B,C,D,E (Art. No.: R4101) is an enzyme immunoassay for the identification of staphylococcal enterotoxins A, B, C, D and E in fluid and solid foods as well as in bacterial cultures.

All reagents required for the enzyme immunoassay are contained in the test kit. The test kit is sufficient for 12 determinations. A microtiter plate spectrophotometer is required for quantification.

Sample preparation: Extraction according to a simplified method for preparation of samples (only conditionally applicable, see section 9. "Sample preparation") by homogenization with buffer and sample preparation

Extraction according to the "Official European Preparation Method" (dialysis concentration method, (see section 9.5)

Time requirement: Sample preparation for 10 samples..... approx. 1 h (simplified preparation)

Sample preparation for 10 samples..... approx. 19 h (official method, dialysis overnight)

Test implementation (incubation time) 2 h 45 min

Limit of detection

Simplified preparation: Liquid samples 0.25 ng/ml toxin
Solid samples..... 0.375 ng/g toxin
Supernatants from bacterial cultures..... 0.25 ng/ml toxin

Dialysis concentration: Liquid samples 0.05 ng/ml toxin
Solid samples..... 0.05 ng/g toxin

Further information is contained in the validation report.

The specificity of the RIDASCREEN® SET A,B,C,D,E test was determined by analyzing the cross-reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice manual. It lists minimum standards concerning the framework conditions when using test kits of R-Biopharm AG and performing ELISA analyses with them. The manual can be retrieved, printed, and downloaded from www.r-biopharm.com/products/food-feed-analysis.

Related products

RIDASCREEN® SET Total (R4105; 96 determinations)

1. Intended use

RIDASCREEN® SET A,B,C,D,E is a sandwich enzyme immunoassay for the identification of staphylococcal enterotoxins (SET) A, B, C, D, E in fluid and solid foods as well as in bacterial cultures. Based on its sensitivity the RIDASCREEN® SET A,B,C,D,E test is consequently clear superior to the immunodiffusion procedure which has a detection limit of 0.1 mg/ml.


2. General information

Staphylococci belong to the family of micrococci. *Staphylococcus aureus* and *Staphylococcus hyicus* produce one or more heat stable proteins that behave as enterotoxins. These are the causative agents for a number of food poisoning cases. The main causative agent of food poisoning are enterotoxins of *Staphylococcus aureus* next to the Salmonella. Generally, it is assumed that a population of 5×10^5 cells of enterotoxin-producing *Staphylococcus aureus* strains per gram of food is required to lead to an intoxication. However, other studies show that only 100 - 200 ng of staphylococcal enterotoxins can lead to symptoms of food poisoning.

SET intoxications have been frequently associated with pasta, finished meat products, ham, pies, chicken meat products, fish, fish products, milk, milk products, ice cream, egg products, salads, pastries and cake fillings as well as preparations from these food products. The enterotoxins of the serological group A, B, C, D, E are very significant.

3. Test principle

The basis of the test is the antigen-antibody reaction. The wells A - E and H in the microtiter strips are coated with specific antibodies against staphylococcal enterotoxins A, B, C, D and E, the wells F and G serve as controls and are coated with antibodies of non-immunized animals in the following pattern:

Microtiter strips of the RIDASCREEN® SET A,B,C,D,E enzyme immunoassay

This sequence is obtained when the microtiter strips are inserted properly into the frame (narrow end of the strip → at the top, wide end of the strip → at the bottom). By adding the sample solution to the wells A to G (see the letters on the left side of the frame) and the positive control to well H, present toxins will bind to specific

capture antibodies. Sample components not bound by the antibodies are then removed in a washing step. The immobilized toxins were bound by a mixture of specific antibodies conjugated to biotin. Any unbound conjugate is then removed in a washing step. After addition of Streptavidin-POD as well as an additional washing step, the detection is performed by adding the substrate/chromogen solution. Bound Streptavidin-POD conjugate converts the reddish chromogen into a blue end product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 / 620 ± 10 nm; the absorbance is proportional to the SET concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for 12 measurements (from the 96 wells of the microtiter plate one strip of 8 single wells is necessary per single test). Each test kit contains:

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		96 wells
Positive control	Red	Ready to use	Red stained	2 ml
Wash buffer	Brown	Concentrate	10 x	100 ml
Conjugate 1	Red	Ready to use	Green stained	11 ml
Conjugate 2	Black	Ready to use	Blue stained	11 ml
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use	Red stained	13 ml
Stop solution	Yellow	Ready to use		14 ml

5. Reagents required but not provided

5.1 Equipment

- Analytical balance and weighing vessels
- Mixer or equivalent for sample homogenization. For food matrices like whey powder and others that are difficult to homogenize, an ultra turrax is highly recommended to obtain an homogeneous sample.
- 50 ml tubes
- Optional: sterile filter
- 100 µl micro pipette
- Incubator 35 - 37 °C (95 - 98.6 °F)
- Multichannel pipette or microwell plate washer
- Microwell plate spectrophotometer (450 / 620 ± 10 nm)

Additional equipment for sample preparation according to the European official screening method (EOSM) Version 5, September 2010:

General remark: it is strictly recommended to use only laboratory vessels (e.g. funnels, beakers, test tubes) of laboratory glass or polypropylene to avoid the absorption of toxins.

- Shaker for beakers (room temperature)
- Refrigerated centrifuge (4 °C), 3130 - 10 000 g, centrifuge tubes
- Dialysis membrane, MWCO: 6 - 8 kD, flat width: 23 ± 2 mm (e.g. Spectra / Por[®]1, ref: 132 650, Spectrum)
- Closures, sealing width = 35 mm (e.g. Spectra / Por[®], ref: 132 736, Spectrum)
- pH-meter
- 50 ml tubes
- Funnel
- Glass wool
- Laboratory glass trough
- Refrigerator (5 ± 3 °C / 41 ± 5.4 °F) and freezer (≤ -18 °C / ≤ -0.4 °F)
- Vortexer
- Microwell plate spectrophotometer (450 / 620 ± 10 nm)

5.2 Reagents

- PBS buffer, pH 7.4, (0.55 g NaH₂PO₄ x H₂O + 2.85 g Na₂HPO₄ x 2 H₂O + 8.7 g NaCl, fill up to 1000 ml with distilled water)
- Distilled water (optional: sterile water)
- n-heptane (for samples with high fat content)
- Brain heart infusion (BHI) for the pre-enrichment of potentially toxin-forming Staphylococci strains. Ask your local producer / distributor of culture media for the supply with BHI.

Additional reagents for EOSM:

- Hydrochloric acid (5N and 1N)
- Sodium hydroxide (5N and 1N)
- Polyethylene glycol 20 000 (PEG), quality for synthesis

6. Warnings and precautions for the users

This test should be carried out only by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. Please refer to the component safety information in the material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

Ensure the proper and responsible disposal of all reagents and materials after their use. For disposal, please adhere to national regulations.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive, therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- Bluish coloration of the reddish substrate/chromogen prior to addition in the wells
- Extinction less than 1.0 ($E_{450/620 \text{ nm}} < 1.0$) for the positive control

9. Sample preparation

The samples have to be stored at 5 ± 3 °C (41 ± 5.4 °F) until extraction. The samples should be completely thawed before the extraction step.

Important notes: To ensure good precipitation of solids and complete phase separation via centrifugation, increase the centrifugation speed and / or time as needed.

Food proteins may interact with toxins or antibodies of the test kit, in which case they can significantly affect the detection of the toxins. Therefore it is strongly recommended to prepare food samples according to the “Official European Screening Method of the European Reference Laboratory for Coagulase-positive Staphylococci” (see 9.5).

The simplified methods for sample preparation described under 9.1 to 9.3 may not be suitable for extracting SET from difficult matrices (e.g. fish, chocolate, pickled vegetables). For these foods, particularly low recovery rates and / or non-specific bindings of matrix proteins to the test antibodies have been observed.

A list of SET recovery rates determined for a wide range of different food matrices is available on request.

9.1 Simplified preparation of milk

- Centrifuge milk samples (10 - 25 ml), especially raw milk samples, in cool condition: 10 min / 3500 g / 10 °C (50 °F)
(if no refrigerated centrifuge is available, pre-cooling of samples is necessary)
- Thoroughly remove the upper cream layer
- Use 100 µl per well in the assay

9.2 Simplified preparation of pasta and rice (cooked), meat, ice cream, processed foods and other foods with a fat content of less than 40 %

- Grind 10 - 25 g of the sample and homogenize with 1.5 ml of PBS buffer (pH 7.4) per g sample (e.g. 10 g sample + 15 ml buffer)
- Shake for 15 min
- Centrifuge: 10 min / 3500 g / 10 °C (50 °F)
- If necessary, remove upper fat layer
- Use 100 µl per well in the assay

9.3 Simplified preparation of foods with a fat content of more than 40 %

- Grind 10 - 25 g of the sample and homogenize with 1.5 ml of PBS buffer (pH 7.4) per g sample (e.g. 10 g sample + 15 ml buffer)
- Shake for 15 min

- Centrifuge: 10 min / 3500 g / 10 °C (50 °F)
- Transfer aqueous phase to another centrifugal vial, add the same volume of n-heptane and mix thoroughly for 5 min
- Centrifuge: 5 min / 3500 g / 10 °C (50 °F)
- Thoroughly remove the upper heptane layer; avoid transferring heptane residues into the wells
- Use 100 µl of the resulting aqueous (lower) phase per well in the assay

9.4 Bacterial cultures

Potentially toxin-forming strains of the species *S. aureus* (as well as *S. hyicus* or *S. intermedius*) have to be pre-cultivated in Brain heart infusion (BHI) to ensure optimal formation of enterotoxins.

Important note: Prior to analysis, please make sure that all strains to be tested are present in pure culture.

- Centrifuge supernatants of microbiological fluid cultures for 5 min / at least 3500 g / 10 °C (50 °F)
- Sterile filtration of the supernatant required because unprecipitated or re-suspended cells may disturb the test reaction
- Use 100 µl of the filtrate per well in the assay

Remark: Background effects which may occur (too high values for the negative controls) can be avoided, if the sterile filtrated supernatant is further diluted with PBS buffer.

This is also valid if OD values measured for detected toxins are outside the linear range (about ≥ 3.0 , please note manufacturer's specifications).

Sample preparations as well as culture supernatants can be stored at 5 ± 3 °C (41 ± 5.4 °F) if the analysis can be performed within 48 h. For longer storage periods the solutions have to be frozen at ≤ -18 °C (≤ -0.4 °F).

Frozen sample preparations or culture supernatants have to be completely defrosted and brought to room temperature before performance of analysis. Implement the once thawed samples into RIDASCREEN®SET A,B,C,D,E immediately. They can't be stored deep frozen again.

Sample preparations or culture supernatants prepared for analysis with RIDASCREEN® SET Total can be used directly for toxin identification with RIDASCREEN® SET A,B,C,D,E.

9.5 Sample preparation according to the European official screening method, Version 5, September 2010

9.5.1 Preparatory steps before toxin extraction

As staphylococcal enterotoxins can be heterogeneously distributed in a sample, mix the whole sample if possible, or a representative part of it, with a mixer.

Weigh $25 \text{ g} \pm 0.1 \text{ g}$ of the homogenized sample and transfer it into a beaker.

Note 1: For cheese with rind, take about 10 % rind and 90 % cheese.

Note 2: For samples in powder form, reconstitute the sample by weighing 12.5 g of sample and 12.5 g of distilled water or follow the manufacturer's instructions (e.g. milk powder as a 10 % solution).

Note 3: In the case of a suspected staphylococcal food poisoning outbreak, a minimal of 12.5 g of food should be tested.

9.5.2 Extraction of enterotoxins

Add 40 ml of warm distilled or deionized water ($38 \text{ °C} \pm 2 \text{ °C}$ / $100.4 \pm 7.2 \text{ °F}$) to the test portion and homogenize the mixture by using a turrax or a blender. After use, rinse the mixing system with distilled water and add the rinse water to the sample being analyzed.

Note 1: For a liquid sample, do not add 40 ml of distilled water.

Shake the sample at room temperature for at least 30 min to homogeneously distribute the toxins.

Acidification step: Acidify the mixture with a few drops of hydrochloric acid in order to obtain a **pH between 3.5 and 4.0**.

Note 2: To prevent denaturation of the enterotoxins during acidification, maintain a **pH between 3.5 and 4.0**. A suitable pH meter must be used for

adjustment. Do not allow pH to fall under 3.0. If the pH drops below 3.0 during acidification, take a fresh 25 g test portion and prepare as described in 9.5.1.

Centrifuge the mixture at least at 3130 x g for 15 min at **4 °C (32.9 °F) or room temperature (20 - 25 °C)**. Transfer the supernatant into a beaker.

Note 3: It is recommended to rinse vessels and devices with distilled water after each step to recover a maximum of toxin.

Note 4: If the supernatant is not clear enough, centrifuge again as described above.

Note 5: The pH of the supernatant after the first centrifugation has to be < 4.5. If this is not the case, acidify to a **pH between 3.5 and 4.0** and centrifuge again.

Neutralization step: Neutralize the mixture with NaOH solution in order to obtain a **pH between 7.4 and 7.6**. Centrifuge again as described. Remove as much of the neutralized aqueous phase (supernatant) as possible and continue to use.

Note 6: To prevent denaturation of the enterotoxins during neutralization maintain a **pH between 7.4 and 7.6**. A suitable pH meter must be used for adjustment.

Do not allow pH to rise above 9.0. If the pH rises above 9.0 during neutralization, take a fresh 25 g test portion and prepare as described in 9.5.1.

9.5.3 Extract concentration via dialysis

For each sample:

- Prepare a 30 % (w/v) PEG solution (30 g PEG / 100 ml distilled water)
- Cut about 50 to 60 cm off a dialysis membrane
- Soak the membrane in distilled water per the manufacturer's instructions (at room temperature for at least 30 minutes)
- Rinse the membrane with distilled water (outside and inside)
- Lock one end of the membrane with a closure, fill it up with the neutralized aqueous phase as prepared in 9.5.2 by using a funnel and a small piece of glass wool to keep out resuspended particles. Lock the other end of the membrane with a second closure.

Note 1: If the sample to analyze contains high amounts of salt or sugar, conduct a dialysis under agitation with 2 L of distilled water two times in one hour.

- Lay the filled dialysis membrane in a laboratory glass trough containing the 30 % (w/v) PEG solution. Allow the extracts to concentrate overnight at 5 ± 3 °C (41 ± 5.4 °F).

Note 2: If the extract is not concentrated enough, increase the dialysis time as needed. It may be necessary to add more PEG powder to the solution.

- Take the dialysis membrane out of the PEG solution and rinse the outside of the membrane with distilled water to remove all traces of PEG.

Remove the concentrated extract using:

- **PBS solution if of extract contains milk or milk products.**
- **Deionized water if extract contains no milk or milk products.**

Rinse the inside of the dialysis membrane well to obtain a final concentrated extract mass between 5.0 g to 5.5 g (maximum 5.8 g for sticky extracts).

When removing the concentrate from the dialysis membrane, it is recommended:

- To rub the insides of the dialysis membrane together in order to loosen the toxins from the membrane walls and to recover the maximum amount of enterotoxins.
- To pour small drops of PBS or deionized water on the inside of the membrane multiple times and rinse repeatedly to recover all enterotoxins present.

Carefully transfer the concentrated extract into a glass vial.

Note 3: If the initial weight of the sample to be tested is lower than 25 g (9.5.1, note 3), the weight ratio of the weighed test portion to the concentrated extract should be 5:1.

For outbreaks of food poisoning or for special studies, the test portion mass may differ from the size of 25 g. The ratio of test portion mass to concentrated extract is shown in the following table:

Test portion mass	Concentrated extract
17.5 g - 25.0 g	3.5 g - 5.0 g (ratio 5:1)
12.5 g - 17.5 g	3.5 g (3.9 g max.)

Note 4: If the concentrated extract is analyzed within 48 h, store it at $5 \pm 3 \text{ }^{\circ}\text{C}$ ($41 \pm 5.4 \text{ }^{\circ}\text{F}$), otherwise store it at $\leq -18 \text{ }^{\circ}\text{C}$ ($\leq -0.4 \text{ }^{\circ}\text{F}$). The extract should be completely thawed and homogenized before testing.

10. Test procedure

10.1 Preliminary comments

Bring all reagents as well as prepared samples to room temperature ($20 - 25 \text{ }^{\circ}\text{C}$ / $68 - 77 \text{ }^{\circ}\text{F}$) before use.

The **wash buffer** is provided as a 10fold concentrate. Before use, the buffer has to be diluted 1:10 (1 + 9) with distilled water (e.g. 10 ml buffer concentrate + 90 ml dist. water). Prior to dilution, dissolve any formed crystals by heating the buffer in a water bath at $37 \text{ }^{\circ}\text{C}$ ($98.6 \text{ }^{\circ}\text{F}$). The diluted and ready-to-use buffer has a shelf life of one week at $2 - 8 \text{ }^{\circ}\text{C}$ ($35 - 46 \text{ }^{\circ}\text{F}$).

Unused reagents should be immediately stored at $2 - 8 \text{ }^{\circ}\text{C}$.

10.2 Test procedure

Carefully follow the recommended washing procedure to obtain unambiguous results. Do not allow wells to dry between work steps.

1. Insert as many microtiter strips into the microwell holder as samples are present. For each sample, all wells of a microtiter strip are required.
2. Add $100 \text{ }\mu\text{l}$ of a sample to wells A to G of one microtiter strip, add $100 \text{ }\mu\text{l}$ of the positive control to well H. Mix gently by shaking the plate manually. Cover the plate and incubate for 1 h at $35 - 37 \text{ }^{\circ}\text{C}$ ($95 - 98.6 \text{ }^{\circ}\text{F}$).
3. Dump the liquid out of the wells into a sink. Tap the microwell holder upside down onto a clean filter towel (three times in a row) to remove all remaining liquid from the wells. Wash each well with $300 \text{ }\mu\text{l}$ of wash buffer (10.1). Empty the wells again and remove all remaining liquid. Repeat the washing step four more times.
4. Pipette $100 \text{ }\mu\text{l}$ of conjugate 1 into each well. Mix gently by shaking the plate manually. Cover the plate and incubate for 1 h at $35 - 37 \text{ }^{\circ}\text{C}$ ($95 - 98.6 \text{ }^{\circ}\text{F}$).
5. Dump the liquid out of the wells into a sink. Tap the microwell holder upside down onto a clean filter towel (three times in a row) to remove all remaining

liquid from the wells. Wash each well with 300 µl of wash buffer (10.1). Empty the wells again and remove all remaining liquid. Repeat the washing step four more times.

6. Pipette 100 µl of conjugate 2 into each well. Mix gently by shaking the plate manually. Cover the plate and incubate for 30 min at 35 - 37 °C (95 - 98.6 °F).
7. Dump the liquid out of the wells into a sink. Tap the microwell holder upside down onto a clean filter towel (three times in a row) to remove all remaining liquid from the wells. Wash each well with 300 µl of wash buffer (10.1). Empty the wells again and remove all remaining liquid. Repeat the washing step four more times.
8. Pipette 100 µl of substrate/chromogen to each well. Mix gently by shaking the plate manually. Cover the plate and incubate for 15 min at 35 - 37 °C (95 - 98.6 °F).
9. Pipette 100 µl of stop solution into each well. Mix gently by shaking the plate manually and measure the extinction at 450 / 620 ± 10 nm in a microtiter plate spectrophotometer. Read within 30 minutes after addition of stop solution.

11. Evaluation

Special software, **RIDASOFT® Win.NET** (Art. No. Z9996FF), is optionally available for evaluating RIDASCREEN® enzyme immunoassays.

11.1 Quality control of the test

- The extinction value for the positive control should be equal to or greater than 1.0.
- The mean extinction value for the negative controls should be equal to or less than 0.2 units.

If these criteria have not been met, the following steps should be checked and corrected before repeating the test:

- Check kit expiry date.
- Ensure sufficient time was allowed for kit components to reach room temperature (20 - 25 °C / 68 - 77 °F).
- A new pipette tip was used for each sample or control to avoid cross-contamination.
- Check sterility of water for preparing wash buffer (use sterile water for preparation).

- Inadequate washing of wells / fewer washing steps than recommended (see 10.2).
- Contaminated pipettes (clean regularly).

11.2 Determination of threshold values

The mean absorbance value for the negative controls in wells F and G is determined separately for each sample. Add 0.15 to the mean value of the negative control to get the threshold value.

Example: Negative control 1 (well F) = 0.008
 Negative control 2 (well G) = 0.010
 Mean value = 0.009
 Threshold value = 0.009 + 0.15 = 0.159

11.3 Result interpretation

- A sample is considered POSITIVE for a certain SET if the test is valid (11.1) and the extinction measured at 450 / 620 ± 10 nm in the appropriate well is greater than or equal to the threshold value.
- A sample is considered NEGATIVE for a certain SET if the test is valid (11.1) and the extinction measured at 450 / 620 ± 10 nm in the appropriate well is less than the threshold value.

11.4 Cross-reactivities of antibodies

Cross-reactivities are known to occur between antibody/toxin: A/E, E/A, B/C and C/B. The percentage of cross-reactivities is between 10 - 50 %. Some food matrices (e.g.: milk, sausages) may enhance or weaken the cross-reactivity. If cross-reactivities are suspected the prepared samples should be diluted with PBS buffer (see below) and measured again.

To distinguish significantly between cross-reactivity and actual detection the absorbance values of positive toxin measurements should be in a range between 1.0 - 2.5. For positive samples with higher extinction values then 2.5 the sample should be diluted, that the specific signal is in the range of 1.0 - 2.5.

Example 1 for **cross-reactivity** between A/E:

A) Initial sample results:

Extinction value for detection of Toxin A: 4.1

Extinction value for detection of Toxin E: 3.6

B) The sample should be diluted to obtain an extinction value in the range between 1.0 and 2.5:

Extinction value for detection of Toxin A: 2.5

Extinction value for detection of Toxin E: **1.0**

Explanation: The signal of the sample in the cavity of toxin E is 1.0 and would be rated as strong positive. But it is only 40 % of the specific signal of the sample in cavity A. This has to be evaluated as cross-reactivity. The sample only contains toxin A.

Example 2 for **no cross-reactivity** between A/E:

A) Initial sample results:

Extinction value for detection of Toxin A: 4.1

Extinction value for detection of Toxin E: 3.6

B) The sample should be diluted to obtain an extinction value in the range between 1.0 and 2.5:

Extinction value for detection of Toxin A: 2.5

Extinction value for detection of Toxin E: **1.9**

Explanation: The main difference between example 1 to example 2 is the extinction value of the sample in cavity E (1.9), which is clearly above 50 % of the specific signal of the sample in cavity A. In conclusion, there is no cross-reactivity. The sample contains toxin A and toxin E.

Further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2017-08-09	Release version
2020-10-15	10.1 Preliminary comments, 11 Evaluation

Explanation of symbols

- General symbols:

	Follow the instructions for use
	Batch number
	Expiry date (YYYY-MM)
	Storage temperature
	Article number
	Number of test determinations
	Manufacturing date (YYYY-MM)
	Manufacturer + address

Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.