# RIDA® STAMP

## Art. No.: HS0291 – HS2012

HS0291; HS0292	RIDA® STAMP Total
HS0371; HS0372	RIDA® STAMP YM-P
HS0392	RIDA® STAMP Salmonella
HS0411; HS0412	RIDA® STAMP Coliform
HS0431; HS0432	RIDA® STAMP ECC
HS0462	RIDA® STAMP S. aureus
HS1831; HS1832	RIDA® STAMP Total Desi
HS2011; HS2012	RIDA® STAMP Pseudomonas

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Ready-to-use agar stamp plate for enumeration of  
microorganisms on surfaces

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Please contact for questions and further information:

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RIDA® and RIDASCREEN®

are registered trademarks of R-Biopharm AG

Manufacturer: R-Biopharm AG, Darmstadt, Germany

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# RIDA<sup>®</sup>STAMP

## 1. Intended use

RIDA<sup>®</sup>STAMPs are ready-to-use agar plates for quantitative detection of microorganisms on surfaces. According to the specification of the individual nutrient medium only the species given in the instructions will be detected. The tests are available in different package sizes of 50 or 25 single tests.

## 2. General

Microorganisms of various species occur in all imaginable habitats, worldwide (ubiquitous organisms). They can be found naturally in soils, on plants, in waters and in the gut of humans and animals and (transferred by dusts and aerosols) are constant parts of the ambient air.

Especially in the process of food production contaminations of surfaces in the production environment are inevitable. Good manufacturing praxis (GMP), well-conceived hygiene management as well as the adherence of general workplace hygiene are crucial factors to decrease the risk of implementing surface contaminations into the food stuff.

Therefore the correct performance of cleaning and sanitation measurements is absolutely important and certainly the regularly verification of successful performed disinfection is important, too.

The microbiological agars of the "RIDA<sup>®</sup>STAMP" product line are perfectly suited for direct stamps of surfaces in the production environment and therefore useable as control tests. RIDA<sup>®</sup>STAMP plates can be applied also on solid foods to make microbial loads on food surfaces visible.

## 3. Test principle

The RIDA<sup>®</sup>STAMP plate consists of a 10 cm<sup>2</sup> petri dish with pedestal. The petri dish contains the individual specific agar according to each parameter. In each agar type specific nutrients, selective agents as well as chromogenic enzyme substrates are implemented depending on its specification. These ingredients are responsible for formation of typical colonies during growth of microorganisms. Subsequently to the defined incubation time the colonies can be counted and the numbers may be taken to rate the degree of contamination on the sampled surfaces.

## 4. Reagents provided

In accordance to the package size agar stamp plates for performance of 50 or 25 determinations are contained. The plates are stretched in bars of 5 and are singly breakable. Each bar is welded in a plastic bag.

An interpretation guide sheet with sample images is enclosed in each package.

Production codes on the plastic wrappings:

RIDA®STAMP Total .....	SMA*
RIDA®STAMP YM-P .....	PDA
RIDA®STAMP Salmonella.....	MCLB
RIDA®STAMP Coliform .....	X-GAL
RIDA®STAMP ECC.....	XM-G
RIDA®STAMP S. aureus .....	XSA
RIDA®STAMP Total Desi .....	SCDLP
RIDA®STAMP Pseudomonas.....	PSEU

\*SMA (Standard Method Agar = Total Plate Count Agar)

## 5. Materials required but not provided

5.1. Equipment (obligatory or optional according to application):

–Incubator (30 - 37 °C / 86 - 98.6 °F)

5.2. Reagents and disposals

For application of RIDA®STAMP direct stamp plates no additional reagents or disposals are necessary.

## 6. Warnings and precautions for the users

Discard the used direct stamp plates according to the current disposal regulations for use of microbiological media (e.g. by autoclaving before handing them over to the waste).

## 7. Storage instructions

Store RIDA®STAMP at 2 - 8 °C (35.6 - 46.4 °F). Do not freeze agar plates in any case!

Return any unused direct stamp plates to their original plastic bag and seal it. Store at 2 - 8 °C (35.6 - 46.4 °F) further on. Use plates from opened bags within a month.

The maximum shelf life of RIDA®STAMP is for:

RIDA®STAMP Total .....	12 months
RIDA®STAMP YM-P .....	12 months
RIDA®STAMP Salmonella.....	5 months
RIDA®STAMP Coliform .....	12 months
RIDA®STAMP ECC.....	12 months
RIDA®STAMP S. aureus .....	4 months
RIDA®STAMP Total Desi .....	12 months
RIDA®STAMP Pseudomonas.....	7 months

## 8. Indication of instability or deterioration of reagents

Changes in coloration of agar (darkening of light-yellowish agars, bleaching of colored agars, etc.) are indicating instability or deterioration of reagents or nutrients. Shrinking of agar area (loss of over-top area of agar, appearing gaps between edge of plate and agar area) indicates water loss/beginning dry out. Do not use discolored or shrunken stamp plates!

## 9. Preparation of samples

The sampling surface should be plane and free of rough particles. Foods are only suitable for direct stamps if their surface is plane, solid and non-adherent. The food surface should not contain any protruding compounds which could damage the surface of the agar (bone fragments, fish bones, etc.).

**Remark:** If disinfectant activity is decreasing, lethally damaged microorganisms may remain on the surfaces. During the direct stamp procedure they get onto the agar together with residues of disinfectants. The disinfectant residues inhibit the growth of the still living organisms and therefore their detection. In case of disinfectant residues or to control the efficiency of disinfectants RIDA®STAMP Total Desi (HS1831/HS1832, contains neutralization agents) should be used.

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- Open plastic bag using scissors
- Take out agar plates and break the number of plates necessary. Put back the remaining plates, seal bag with tape and store at 2 - 8 °C (35.6 - 46.4 °F)
- Mark stamp plate with identity of sampling sites

- Open lid of stamp plate and press plate against sampling side with gentle pressure
- Close stamp plate and incubate at the required temperature (10.2.)

## 10.2. Incubation

Turn the direct stamp plates upside down and stack them sorted according to plate types (there might be different incubation temperatures for each plate type, s. below). RIDA<sup>®</sup>STAMP direct stamp plates are incubated under the temperature conditions mentioned in the following:

RIDA <sup>®</sup> STAMP Total .....	30 - 35 °C (86 - 95 °F) / 24 - 48 h
RIDA <sup>®</sup> STAMP YM-P .....	20 - 25 °C (68 – 77 °F) / 2 – 5 d
RIDA <sup>®</sup> STAMP Salmonella.....	36 ± 1 °C (96.8 ± 1.8 °F) / 24 – 48 h
RIDA <sup>®</sup> STAMP Coliform .....	36 ± 1 °C (96.8 ± 1.8 °F) / 24 – 48 h
RIDA <sup>®</sup> STAMP ECC.....	35 ± 0.5 °C (95 ± 0.9 °F) / 20 ± 2 h
RIDA <sup>®</sup> STAMP S. aureus .....	36 ± 1 °C (96.8 ± 1.8 °F) / 22 – 24 h
RIDA <sup>®</sup> STAMP Total Desi .....	30 - 35 °C (86 - 95 °F) / 24 - 48 h
RIDA <sup>®</sup> STAMP Pseudomonas.....	35 ± 2°C (95 ± 3.6 °F) / 24 - 48 h

## 11. Results

Count the specific colonies which are grown on the RIDA<sup>®</sup>STAMP agar surface after incubation. Classify the results according to your specifications. The interpretation guide sheet contains sample images to be used to compare with colonies on direct stamps.

### 11.1. Individual evaluation of RIDA<sup>®</sup>STAMP Total

Growing microorganisms will mostly form milky-white colonies on the agar surface. Organisms can also be able to produce yellow, orange or red pigments or might be completely colorless. Therefore colored or transparent colonies may appear as well.

Count all visible colonies on the agar surface.

### 11.2. Individual evaluation of RIDA<sup>®</sup>STAMP YM-P

Yeasts form almost small round shaped to oval colonies which are normally of whitish to rose color. Molds form their typical air mycelia which can be white brown, green grey or black depending on the individual genus of fungi. All colonies on the agar surface have to be counted to determine the total yeast & mold count.

### 11.3. Individual evaluation of RIDA®STAMP Salmonella

*Salmonella* are producing hydrogen sulfide during growth. In case of presence of iron-(III)-ions precipitates of black iron sulfide are formed. These will turn the whole colonies black or build a black center in the colonies. Due to acid production the agar may change its color from purple to yellow beige.

Different species of *Citrobacter* may form black colonies, too and can be difficult to differentiate.

Reddish or otherwise colored colonies might be formed by other enterobacteria, but are never any kind of indication of *Salmonella* presence.

### 11.4. Individual evaluation of RIDA®STAMP Coliform

Coliform bacteria have the enzyme  $\beta$ -galactosidase which is able to cleave the chromogenic substrate X-GAL contained in the agar. Therefore coliforms form typical blue/blue-green colonies. Other bacterial species are not able to grow or their colonies are of different coloration.

All grown blue or blue-green colonies have to be counted to determine the number of total coliform bacteria.

### 11.5. Individual evaluation of RIDA®STAMP ECC

Due to the combination of the chromogenic substrates X-GLUC and MAGENTA-GAL in the agar, *E. coli* develops blue to blue-red colonies while all other coliform bacteria form red-violet or pink colored colonies.

Nearly all strains of the *E. coli* serovar O:157 do not have any  $\beta$ -Glucuronidase and are therefore not able to form blue colonies. This is why they appear in red-violet and are not to differentiate from other coliforms.

Extended incubation time may promote the growth of accompanying bacteria. If lactic acid bacteria are developing on the agar, these may also cause red-violet or pink colored colonies (the same looking as coliforms), due to their own  $\beta$ -galactosidase activity.

All red-violet or pink colored colonies as well as all blue/blue-red colonies have to be counted to determine the total number of coliform bacteria.

### 11.6. Individual evaluation of RIDA®STAMP S. aureus

*Staphylococcus aureus* forms distinct blue, concave colonies on the agar surface. Coagulase-negative Staphylococci form small colonies appearing in white or blue. Bacteria of the genus *Bacillus* eventually produce faint blue colonies which are flat and lusterless and are clearly to distinguish from real *S. aureus* colonies.

All colonies specified as *S. aureus* have to be counted to determine the total amount of potentially pathogenic Staphylococci on the sampled surface.

#### 11.7. Individual evaluation of RIDA®STAMP Total Desi

According to the description given for the total plate count agar, RIDA®STAMP Total, all colonies grown on the surface of Total Desi have to be counted in the same way.

#### 11.8. Individual evaluation of RIDA®STAMP Pseudomonas

Most *Pseudomonas* species (inclusively *P. aeruginosa*) develop round compressed or dome-shaped colonies with yellow-greenish fluorescence on the medium.

Count all appearing colonies due to the given specification.

### 12. Further analysis

Single colonies can be taken from the agar surface using a loop and can be transferred to other agar or liquid media suitable for subculture or confirmation. Please mind that microbiological work with potential pathogenic organisms (e.g. transfer of *Salmonella* to confirmation media) has to be done only in a suitable microbiological laboratory by trained professionals.

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.