**CONGEN**

# **SureFast® STEC 4plex ONE**

Art. No. F5265  
100 rxn

## **User Manual**

**September 2022**

 **Inhalt**

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## **Content**

1	General Information .....	10
1.1	Description .....	10
1.2	Limit of Detection .....	10
1.3	Kit components and storage .....	11
1.4	Additionally required equipment and materials .....	11
1.5	Setup .....	11
1.6	Detection channel Set-up .....	12
2	Qualitative Analysis .....	13
2.1	Protocol .....	13
2.1.1	Cultural Enrichment .....	13
2.1.2	DNA-preparation .....	13
2.1.3	Preparation of the master-mix .....	14
2.1.4	Preparation of the real-time PCR-mix .....	14
2.2	Interpretation of results .....	15
3	Further Information .....	16
3.1	Product Information .....	16
3.2	Technical Support .....	16

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## **1 General Information**

### **1.1 Description**

The SureFast® STEC 4plex ONE can be applied for the fast and simple isolation, detection and differentiation of specific DNA sequences of the *Escherichia coli* virulence factors *stx1* (subtype a-d), *stx2* (subtype a-g), *eae* and the *Escherichia coli* serotype O157 in different types of food.

The method consists of the following three steps: 1) cultural enrichment, 2) DNA-preparation, 3) specific real-time PCR detection and 4) interpretation of results.

Each reaction contains an internal amplification control (IAC). If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). In addition spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of four fluorescence emissions at the channels FAM, VIC/HEX, ROX und Cy5 at the same time. The technical verification of instruments was performed on Agilent Mx3005P, Agilent AriaDx, Bio-Rad CFX96, Roche LightCycler® 480 II, Roche cobas® z 480 Analyzer, Applied Biosystems 7500, Qiagen Rotor-Gene Q and Bio Molecular Systems MIC.

### **1.2 Limit of Detection**

The SureFast® STEC 4plex ONE real-time PCR has a limit of detection of  $\leq 5$  DNA copies.

The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

The SureFast® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

**1.3 Kit components and storage**

<b>Kit Code</b>	<b>Reagent</b>	<b>Amount</b>	<b>Lid Color</b>
1	Reaction Mix	2 x 1100 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue
L	Lysis Buffer	2 x 25 ml	Clear

**Store all reagents at –20°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.**

The Lysis Buffer should be stored at -20 to +8 °C. Reagents should be used until the expiry date indicated on the label.

**Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.**

**1.4 Additionally required equipment and materials**

- real-time PCR instrument with four detection channels (510 nm, 580 nm, 610 nm and 660 nm)
- real-time PCR consumable (plates, tubes, capillaries, foils, caps)
- pipettes with filter tips
- 2 ml tubes
- powder-free disposable gloves
- Vortex mixer
- centrifuge with a rotor for the reaction tubes
- Thermomixer/heating block (up to 95 °C)
- enrichment broth (e. g. modified tryptone soy broth (mTSB) with Novobiocin)

**1.5 Setup**

	<b>Blockcycler &amp; R-Biopharm RIDA®CYCLER</b>	<b>Rotorcycler &amp; LightCycler® 480 II</b>
Initial Denaturation (HOLD) Cycles	1 min, 95°C 45	1 min, 95°C 45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

**1.6 Detection channel Set-up**

Real-time PCR device	Detection	Detection channel	Quencher	Note
<b>Agilent Mx3005P</b>	<i>E. coli stx1/stx2</i>	FAM	+	
	IAC	HEX	+	
	<i>E. coli</i> O157	ROX	+	
	<i>E. coli eae</i>	Cy5	+	
<b>Agilent AriaDx</b>	<i>E. coli stx1/stx2</i>	FAM	+	
	IAC	HEX	+	
	<i>E. coli</i> O157	ROX	+	
	<i>E. coli eae</i>	Cy5	+	
<b>Applied Biosystems 7500</b>	<i>E. coli stx1/stx2</i>	FAM	None	Check the passive reference option ROX is none.
	IAC	HEX	None	
	<i>E. coli</i> O157	ROX	None	
	<i>E. coli eae</i>	Cy5	None	
<b>Bio-Rad CFX96</b>	<i>E. coli stx1/stx2</i>	FAM	+	
	IAC	VIC/HEX	+	
	<i>E. coli</i> O157	ROX	+	
	<i>E. coli eae</i>	Cy5	+	
<b>R-Biopharm RIDA®CYCLER</b>	<i>E. coli stx1/stx2</i>	green	+	
	IAC	yellow	+	
	<i>E. coli</i> O157	orange	+	
	<i>E. coli eae</i>	red	+	
<b>Qiagen Rotor-Gene Q</b>	<i>E. coli stx1/stx2</i>	green	+	<b>Note:</b> Please use only 0.1 ml reaction tubes. The gain settings must be set to 5 (factory default) for all channels.
	IAC	yellow	+	
	<i>E. coli</i> O157	orange	+	
	<i>E. coli eae</i>	red	+	
<b>Roche LightCycler® 480 II</b>	<i>E. coli stx1/stx2</i>	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	IAC	533-580	+	
	<i>E. coli</i> O157	533-610	+	
	<i>E. coli eae</i>	618-660	+	
<b>Roche cobas® z 480 Analyzer</b>	<i>E. coli stx1/stx2</i>	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	IAC	540-580	+	
	<i>E. coli</i> O157	540-610	+	
	<i>E. coli eae</i>	610-670	+	

## **2 Qualitative Analysis**

### **2.1 Protocol**

#### **2.1.1 Cultural Enrichment**

Sample preparation of food and feed matrices should apply to ISO 6887 and ISO 13136 standards. For enrichment, a volume of 25 g or 25 ml sample material is homogenized with 225 ml of an enrichment broth (e. g. modified tryptone soy broth (mTSB) with Novobiocin). If test portions deviate from 25 g, the enrichment broth volume should be chosen in a 1:10 ratio. Enrichment is performed for 18 h - 24 h at 37 °C ± 1°C. Alternatively other suitable, validated enrichment procedures can be used.

To assess the process of bacterial growth, it is recommended to compare samples at the beginning and at the end of the culturing (bacterial growth at Cp difference > 3).

#### **2.1.2 DNA-preparation**

If the enrichment culture has been shaken or stirred, let the solids settle out for 5 to 10 min.

Transfer 200 µl of the upper third from the enrichment culture into a 2.0 ml tube (not provided with the kit) and add 500 µl of Lysis Buffer (**Code L**) to this sample.

Vortex the tube briefly and put the tube in the heating block.

Incubate for 10 minutes at 95 °C without shaking.

Remove the tube from the heating block and leave the tube at room temperature for 1 minute.

**Note: It is necessary to ensure that the sediment is not stirred up and no particles are pipetted in the PCR reaction. If the particles do not sediment, the lysate should be centrifuged for 1 minute at 12.000 rpm. Afterwards transfer 100 µl of the supernatant into a new tube (not provided with the kit).**

The lysate is ready-to-use for PCR. If the lysate is not used within four hours in the PCR or is intended for a longer storage, transfer 100 µl of the lysate in a new tube (not provided with the kit) and store at -20 °C.

**2.1.3 Preparation of the master-mix**

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay as well as for the inhibition control.

Recommended control reactions for the specific PCR assay: negative control, extraction control, positive control. For the analysis of enrichments additional controls are recommended: zero control (sample before enrichment) and medium control. The reaction mix contains an internal amplification control (IAC) per reaction.

**Reactions needed for the qualitative STEC detection:**

3 reactions for controls (1x no-template control, 1x extraction control, 1x positive control)

For each sample: at least 1 reaction for each sample DNA

**Reactions needed for the qualitative STEC detection in enrichments:**

5 reactions for controls (1x no-template control, 1x extraction control, 1x positive control, 1x zero control, 1x medium control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

**Example for the calculation and preparation of 10 reactions:**

<b>Components of the master-mix</b>	<b>Amount per reaction</b>	<b>10 reactions (with 10% excess)</b>
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

**Mix each master-mix well and centrifuge shortly before use.**

**2.1.4 Preparation of the real-time PCR-mix**

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

**2.2 Interpretation of results**

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

*Escherichia coli stx1/stx2* DNA is detected in the FAM-channel, *Escherichia coli* O157 DNA is detected in the ROX-channel and *Escherichia coli eae* DNA is detected in the Cy5-channel (see table). In the VIC/HEX-channel the amplification control is detected.

A sample is stated **positive** for the respective parameter, if the sample DNA shows amplification in the respective channel. High amplicon concentrations can result in a weak or absent signal of the internal amplification control (IAC).

A Cp value for the internal amplification control (IAC) is not needed to obtain a positive result of the positive control.

A sample is stated **negative** for the respective parameter, if the sample DNA shows no amplification in the respective channel and if the internal control (VIC/HEX-channel) of the sample is **positive** with a shift in Cp-Value  $\leq 2$  compared to the negative control. If the sample DNA in the VIC/HEX channel shows **no amplification** or a shift in Cp-value  $> 2$  compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific *Escherichia coli stx1/stx2*, *Escherichia coli* O157 or *Escherichia coli eae* PCR assay.

Result in the respective channel				Interpretation
FAM channel <i>Escherichia coli stx1/stx2</i>	ROX channel <i>Escherichia coli</i> O157	Cy5 channel <i>Escherichia coli eae</i>	VIC/HEX channel IAC	
<b>positive</b>	negative	negative	<b>positive/negative</b>	<i>Escherichia coli stx1/stx2</i> DNA detected
negative	<b>positive</b>	negative	<b>positive/negative</b>	<i>Escherichia coli</i> O157 DNA detected
negative	negative	<b>positive</b>	<b>positive/negative</b>	<i>Escherichia coli eae</i> DNA detected
negative	negative	negative	positive	Negative, <i>Escherichia coli stx1/stx2</i> , O157, <i>eae</i> – DNA is not detected
negative	negative	negative	negative	invalid

### **3 Further Information**

#### **3.1 Product Information**

- Detailed information about setup of several real-time PCR devices  
(Download: [www.congen.de/en/company/downloads](http://www.congen.de/en/company/downloads))
- Validation Report upon request

#### **3.2 Technical Support**

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de).